

Closed IV System Product Listing

Closed IV system are designed to be used when delivering TPN (Total Parenteral Nutrition) and other medications. The intention of these sets is to reduce the number of touch contaminations that can occur in the NICU. By assembling all pieces together and adding the unique combination of stopcock and needleless valves to our IV medication tubing, we have created an improved delivery system for your IV fluid administration.


BENEFITS (depending on specific set and hospital protocol):

- Reduces the number of connections NICU nurses need to make when setting up IV to deliver TPN and medications, there by reducing chance for infection.
- Stopcock use in combination with a larger syringe of saline will reduce the number of times the system is accessed by 50% or more.
- Needleless valve keeps system closed even when the medication syringe is disconnected.
- Use of color coded clamps reduce the chance of mixing up an IV line.
- Medsafe™ connectors reduce priming volumes as well as minimizing the exposure of incompatible medications.
- Saves money by reducing the number of IV flush syringes used.

CATALOG NUMBER	DESCRIPTION	QTY/BOX
----------------	-------------	---------


AG-089-01	Minibore tubing, vented spike set, back check valve, 4-way stopcock, 0.2 micron filter, tri-furcated extension set.	25
-----------	---	----

Priming Volume: 4.8mL, Length: 79"




AG-100-01	Low volume minibore tubing, vented spike set, pinch clamp, back check valve, 4-way stopcock, male luer lock	50
-----------	---	----

Priming Volume: 7.5ml, Length: 100"



EB-062-01	Microbore tubing, back check valve, 4-way stopcock, female luer locks, slide clamps, 0.2 micron filter, rotating male luer lock	50
-----------	---	----

Priming Volume: 1.4ml Length: 62"



Acacia™


Ph 714.257.0470 Fx 714.257.0513 • www.acaciainc.com

Closed IV System Product Listing

CATALOG NUMBER	DESCRIPTION	QTY/BOX
----------------	-------------	---------


HR-004	Microbore tubing, 4-way largebore stopcock, female luer locks, non-removable slide clamps, 0.2 micron filter, rotating male luer lock, Needleless valve.	50
--------	--	----

Priming Volume: 1.4ml Length: 62"



HR-006	Macrobores tubing, large vented spike, pinch clamp, back check valve, 4-way largebore stopcock, needleless valve, microbore tubing, female luer locks, non-removable slide clamps, 0.2 micron filter, mini tri connector with anti-siphon valve	50
--------	---	----

Length: 95"




HR-008	Microbore tubing, 4-way largebore stopcock, non-removable slide clamp, rotating male luer lock, Needleless valves.	50
--------	--	----

Priming Volume: 0.65ml Length: 62"




HR-009	Microbore tubing, 4-way largebore stopcock, needleless valves, non-removable slide clamps, 0.2 micron filter, male luer lock, back check valves, rotating male luer lock	50
--------	--	----

Priming Volume: 1.75ml Length: 62"



HR-012	Microbore tubing, large vented spike, pinch clamp, check valves, needleless valve, non-removable slide clamp, anti-siphon valve. (closed system)	50
--------	--	----

Length: 95"



Acacia™

Ph 714.257.0470 Fx 714.257.0513 • www.acaciainc.com

IMA
Independent Medical Associates

Products distributed by Independent Medical Associates: www.I-MA.com · info@I-MA.com · p: 888.548.4462 · f: 888-548-1462