

INTOPO™ Enteral Feeding System

How Safe are Your Enteral Feeding Lines?

Tubing misconnections can have devastating effects for the patient. Unfortunately, these misconnections do happen. Since caring for your patient is priority number one, shouldn't you also be taking measures to ensure your patient's safety?

Introducing the INTOPO™ Enteral Feeding System

Our enteral products are designed to prevent potentially deadly tubing misconnections. With color coded connectors and tubing, you'll be able to clearly and immediately identify enteral feeding lines. Based on recommendations by JCAHO's Sentinel Alert, our connectors and tubing will not physically mate with IV luers.

A COMPLETE SYSTEM: Don't allow yourself the false security of an incomplete system. Our orange color coded connectors and tubing start at the feeding line and end with the syringe.

PUMPING ACCURACY: Our syringes are compatible with pumps so you'll know your patient is getting the right amount every time.

Safety Features and Benefits

- Color Coded: Orange for clear and immediate visual distinction and identification
- Compatible: Syringes are compatible and accurate with pump use
- Connectability: Will not connect to any standard luer



Set change recommendation is every 4 hours.

"Hang times for formula in the NICU, for other immune-suppressed patients, and for human milk should be 4 hours or less." From the *Guidelines for the Preparation of Formula and Breast Milk in Health Care Facilities* Pediatric Nutrition Practice Group of the American Dietetic Association, 2003. To view guidelines visit: http://www.eatright.org/Public/NutritionInformation/92_17241.cfm

INTOPO™ Product Listing



INTOPO Syringe: compatible with pumps

CATALOG NUMBER	DESCRIPTION	QTY/BOX
INT-03	3cc syringe	100
INT-05	5cc syringe	100
INT-10	10cc syringe	100
INT-20	20cc syringe	100
INT-30	30cc syringe	50
INT-60	60cc syringe	50





INTOPO Feeding Tube

CATALOG NUMBER	DESCRIPTION	QTY/BOX
ENPU-20-50	5 French, 20" Polyurethane feeding tube for use up to 30 days	10
ENPU-36-50	5 French, 36" Polyurethane feeding tube for use up to 30 days	10
ENPU-20-80	8 French, 20" Polyurethane feeding tube for use up to 30 days	10
ENPV-16-50	5 French, 16" PVC feeding tube for use up to 3 days	50
ENPV-16-80	8 French, 16" PVC feeding tube for use up to 3 days	50

EXTENSION SETS

CATALOG NUMBER	DESCRIPTION	QTY/BOX
EB-060-18	Orange-striped low volume minibore tubing, oral female adapter, non-removable slide clamp, universal enteral connector. Priming Volume: 1.3mL, Length: 60"	100
EB-060-19	Orange-striped low volume minibore tubing, oral female adapters, non-removable slide clamps, non-vented tethered cap, universal enteral connector. Priming Volume: 1.7mL, Length: 60"	100



Acacia™

Ph 714.257.0470 Fx 714.257.0513 • www.acaciainc.com