

# INTOPO™ Enteral Feeding System

## How Safe are Your Enteral Feeding Lines?

Tubing misconnections can have devastating effects for the patient. Unfortunately, these misconnections do happen. Since caring for your patient is priority number one, shouldn't you also be taking measures to ensure your patient's safety?



## Introducing the INTOPO™ Enteral Feeding System

Our enteral products are designed to prevent potentially deadly tubing misconnections. With color coded connectors and tubing, you'll be able to clearly and immediately identify enteral feeding lines. Based on recommendations by JCAHO's Sentinel Alert, our connectors and tubing will not physically mate with IV luers.

**A COMPLETE SYSTEM:** Don't allow yourself the false security of an incomplete system. Our orange color coded connectors and tubing start at the feeding line and end with the syringe.

**PUMPING ACCURACY:** Our syringes are compatible with pumps so you'll know your patient is getting the right amount every time.

## Safety Features and Benefits

- **Color Coded:** Orange for clear and immediate visual distinction and identification
- **Compatible:** Syringes are compatible and accurate with pump use
- **Connectability:** Will not connect to any standard luer



Acacia™

Ph 714.257.0470 Fx 714.257.0513 • [www.acaciainc.com](http://www.acaciainc.com)

**IMA**  
Independent Medical Associates

Products distributed by Independent Medical Associates: [www.I-MA.com](http://www.I-MA.com) · [info@I-MA.com](mailto:info@I-MA.com) · p: 888.548.4462 · f: 888-548-1462