


# MedSafeDemo

## Standard Multi-Access Extension Set

Quadfurcated extension set, minibore tubing, female luer locks, non-removable slide clamps, rotating male luer lock


Approx. Priming Volume: Varies    Approx. Length: 7"



## Acacia's MedSafe Extension Set #BG-005-01

Quadfurcated extension set, minibore tubing, check valves, female luer lock, non-removable slide clamps, rotating male luer lock

Approx. Priming Volume: 0.65ml    Approx. Length: 5"



Acacia®

Ph 714.257.0470    Fx 714.257.0513    [www.acaciainc.com](http://www.acaciainc.com)