

MEDSAFE AND MULTI-ACCESS SETS



BG-004-03
MedSafe Quadfurcated
Extension Set

IMPROVE MEDICATION SAFETY AT THE POINT OF CARE

Acacia Neonatal's MedSafe and Multi-Access Sets are designed specifically to meet the challenges of drug delivery in neonates. With a limited number of viable vascular access points, a multi-route set is essential. Too often, the need arises to create multi-access sets by building awkward and heavy piecemeal sets at the bedside. This practice increases the chance of infection through touch contamination as well as adding to


variability in line set-up, which increases risk of error. Also, the additional tubing length of multiple sets makes delivery of the drug inaccurate and may increase incompatibility issues.

MedSafe and Multi-Access Sets are pre-assembled and patient ready. The MedSafe's unique feature minimizes drug interaction to decrease the chance of dangerous medication errors.

MEDSAFE AND MULTI-ACCESS SETS

PRODUCT FEATURES & BENEFITS:

- **Minimal Dead Space** – The MedSafe connector minimizes potential drug incompatibility in the fluid path, with the shortest common fluid juncture available today.
- **Proper Drug Delivery For Neonates** – MedSafe and Multi-Access Sets are designed for multiple small volume medication delivery, and minimize the risk for incompatibility problems and pooling of medication in the line while minimizing the need for additional flush volume. These features ensure the patient is receiving the medication as administered.
- **Reduced Chance of Infections** – These sets are pre-assembled and patient ready, minimizing the need for additional sets, thereby reducing the risk of infection by touch contamination.
- **Low Priming Volumes** – MedSafe and Multi-Access Sets have very low priming volumes and short lengths.
- **Lightweight** – The specially designed connector is nearly weightless with an easy rotating luer lock to reduce the chance for kinking and pulling.
- **Color Coded** – Slide clamps are different colors for easy, at a glance, line identification.
- **Cost Effective** – Purchasing pre-assembled MedSafe and Multi-Access Sets reduce the need to stock several individual sets leading to a possible decrease in costs.
- **Non-DEHP, Latex Free.**



The MedSafe connector minimizes the sets' dead space.



The sets are pre-assembled and patient ready in Bi-, Tri- (pictured) and quadfurcated configurations.




Each set features color coded slide clamps for ease of use.

BB-004-01

MedSafe Multi-Access Set
Bifurcated
Minibore Tubing
Rotating Male Luer Lock
Two Non-Removable Slide Clamps (White)
Back Check Valve
Female Luer Lock

APPROX. PRIMING VOLUME:	0.4ml
APPROX. LENGTH:	4.0"

50
BOX




BG-004-02

MedSafe Multi-Access Set
Trifurcated
Low-Volume Minibore Tubing
Rotating Male Luer Lock
Three Non-Removable Slide Clamps
(White, Yellow, Blue)
Three Female Luer Locks

APPROX. PRIMING VOLUME:	0.4ml
APPROX. LENGTH:	4.0"

50
BOX




BG-004-03

MedSafe Multi-Access Set
Quadfurcated
Low-Volume Minibore Tubing
Rotating Male Luer Lock
Four Non-Removable Slide Clamps
(White, Yellow, Green, Blue)
Four Female Luer Locks

APPROX. PRIMING VOLUME:	0.5ml
APPROX. LENGTH:	4.0"

50
BOX




BG-004-05

MedSafe Multi-Access Set
Trifurcated
Low-Volume Minibore Tubing
Rotating Male Luer Lock
Three Non-Removable Slide Clamps
(Red, White, Yellow)
Female Luer Lock
Two Back Check Valves

APPROX. PRIMING VOLUME:	0.55ml
APPROX. LENGTH:	4.0"

50
BOX




BG-005-01

MedSafe Multi-Access Set
Quadfurcated
Low-Volume Minibore Tubing
Four Non-Removable Slide Clamps
(White, Yellow, Green, Blue)
Three Back Check Valves
Female Luer Lock

APPROX. PRIMING VOLUME:	0.65ml
APPROX. LENGTH:	5.0

50
BOX




BG-008-01

MedSafe Multi-Access Set
Bifurcated
Minibore Tubing
Rotating Male Luer Lock
0.2 Micron Filter (24 Hours)
Two Non-Removable Slide Clamps
(White, Yellow)
Two Female Luer Locks

APPROX. PRIMING VOLUME:	0.8ml
APPROX. LENGTH:	8.0"

50
BOX




EG-016-01

MedSafe Multi-Access Set
Trifurcated
Low-Volume Minibore Tubing
Rotating Male Luer Lock
Three Non-Removable Slide Clamps
(Yellow, White, Blue)
0.2 Micron Filter (24 Hours)
Three Female Luer Locks
Two Non-Vented Caps

APPROX. PRIMING VOLUME:	0.9ml
APPROX. LENGTH:	16.0"

50
BOX



Distributed By:



11733 66th St N Suite 113
Largo, FL 33773

PH: 888.548.4462
FX: 888.548.1462

E-Mail: info@I-MA.com
Web: www.I-MA.com

© 2011 Acacia, Inc.
Acacia is a registered trademark of Acacia, Inc.

M-0097
021811