













NICU Product Listing

ADMINISTRATION/EXTENSION SETS

CATALOG NUMBER	DESCRIPTION	QTY/BOX	
AG-089-01	Minibore tubing, vented spike transfer set, back check valve, 4-way stopcock, 0.2 micron filter, bifurcated extension set, back check valve. Priming Volume: 4.80mL, Length: 79"	25	
AG-100-01	Low volume minibore tubing, vented spike transfer set, pinch clamp, back check valve, 4-way stopcock, male luer lock. Priming Volume: 7.5mL, Length: 100"	50	
BB-004-01	Minibore tubing, back check valve, female luer lock, two non-removable slide clamps, mini bi-connector w/rotating lock. (Med-Safe) Priming Volume: 0.4mL, Length: 4"	50	
BB-004-03	Med-Safe™, bifurcated extension set, female luer locks, non-removable slide clamps, rotating male luer lock. Priming Volume: 0.4mL, Length: 4"	50	
BB-005-01	Med-Safe™, minibore tubing, bifurcated set, non-removable slide clamps, rotating male luer lock, two needleless valves. Priming Volume: 0.6mL, Length: 5"	50	
BB-006-01	Minibore tubing, female luer locks, removable slide clamps, mini bi-connector, rotating male luer lock. Priming Volume: 0.3mL, Length: 6"	50	







NICU Product Listing

EXTENSION SETS

CATALOG NUMBER	DESCRIPTION	QTY/BOX	
BC-015-02	Kink resistant tubing, bifurcated extension set, pinch clamps, back check valves, rotating male luer lock. Priming Volume: 1.5mL, Length: 16"	50	
BG-004-02	Med-Safe™, low volume minibore tubing, trifurcated extension set, female luer locks, non-removable slide clamps, rotating male luer lock. Priming Volume: 0.4mL, Length: 4"	50	
BG-004-03	Med-Safe™, low volume minibore tubing, female luer locks, non-removable slide clamps, rotating male luer lock. Priming Volume: 0.5mL, Length: 4"	50	
BG-004-04	Med-Safe™, minibore tubing, quadfurcated set, female luer lock, non-removable slide clamps, rotating male luer lock, three needleless valves. Priming Volume: 0.8mL, Length: 5"	50	
BG-004-05	Med-Safe™, minibore tubing, trifurcated extension set, female luer lock, two check valves, non-removable slide clamps, fixed male luer lock. Priming Volume: 0.55mL, Length: 4"	50	
BG-005-01	Med-Safe™, quadfurcated extension set, minibore tubing, check valves, female luer lock, non-removable slide clamps, rotating male luer lock. Priming Volume: 0.65mL, Length: 5"	50	







NICU Product Listing

EXTENSION SETS

CATALOG NUMBER	DESCRIPTION	QTY/BOX	
BG-005-02	Med-Safe™, trifurcated extension set, minibore tubing, check valves, female luer lock, non-removable slide clamps, rotating male luer lock. Priming Volume: 0.5mL, Length: 5"	50	
BG-006-01	Low volume minibore tubing, trifurcated extension set, female luer locks, removable slide clamps, rotating male luer lock. Priming Volume: 0.5mL, Length: 7"	50	
BG-006-02	Low volume minibore tubing, quad-furcated extension set, female luer locks, removable slide clamps, rotating male luer lock. Priming Volume: 0.5mL, Length: 7"	50	
DB-006-17	Med-Safe™, minibore tubing, check valves, mini bi-connector w/ rotating male luer lock Priming Volume: 0.5mL, Length: 6"	50	
DB-010-01	Med-Safe™, minibore tubing, trifurcated set, 0.2 micron filter, female luer locks, non-removable slide clamps, rotating male luer lock, Posiflow™ Needleless Valve. Priming Volume: 3.4mL, Length: 10"	50	
EA-060-01	Microbore tubing, female luer lock, removable slide clamp, male luer lock. Priming Volume: 0.3mL, Length: 60"	50	

NICU Product Listing

EXTENSION SETS

CATALOG NUMBER	DESCRIPTION	QTY/BOX	
EG-016-01	Med-Safe™, low volume minibore tubing, trifurcated extension set, female luer locks, 0.2 micron filter, non-removable slide clamps, rotating male luer lock. Priming Volume: 0.9mL, Length: 16"	50	
EG-016-02	Low-Volume Minibore tubing, female luer lock, 0.2 micron filter, non-removable slide clamps, three check valves, rotating male luer lock Priming Volume: 1ml, Length: 16"	50	
EI-060-02	Large microbore tubing, female luer lock, non-removable slide clamp, male luer lock. Priming Volume: 0.6mL, Length: 60"	50	
PF-019	Med-Safe™, low volume minibore tubing, trifurcated set, female luer locks, 0.2 micron filter, non-removable slide clamps, rotating male luer lock, Posiflow™ Needleless Valve. Priming Volume: 3.6mL, Length: 16"	50	
PF-021	Minibore tubing, non-removable slide clamp, rotating male luer lock, Posiflow™ Needleless Valve. Priming Volume: 0.15mL, Length: 4"	50	
PF-022	Med-Safe™, minibore tubing, quadfurcated set, female luer locks, non-removable slide clamps, rotating male luer lock, Posiflow™ Needleless Valve. Priming Volume: 0.6mL, Length: 4"	50	


NICU Product Listing

EXTENSION SETS

CATALOG NUMBER	DESCRIPTION	QTY/BOX
----------------	-------------	---------


PF-023	Med-Safe™, low volume minibore tubing, trifurcated set, female luer locks, non-removable slide clamps, rotating male luer lock, Posiflow™ Needleless Valve.	50
--------	---	----

Priming Volume: 0.5mL, Length: 4"




PF-024	Med-Safe™, minibore tubing, trifurcated set, female luer lock, 0.2 micron filter, two check valves, non-removable slide clamps, rotating male luer lock, Posiflow™ Needleless Valve.	50
--------	--	----

Priming Volume: 4.6mL, Length: 17"




PF-025	Med-Safe™, minibore tubing, bifurcated set, female luer lock, non-removable slide clamps, rotating male luer lock, Posiflow™ Needleless Valve.	50
--------	--	----

Priming Volume: 0.5mL, Length: 5"




PF-026	Med-Safe™, low volume minibore tubing, trifurcated set, female luer locks, non-removable slide clamps, rotating male luer lock, three Posiflow™ Needleless Valves	50
--------	---	----

Priming Volume: 0.7mL, Length: 5"




PF-027	Med-Safe™, minibore tubing, trifurcated set, female luer lock, 0.2 micron filter, check valve, non-removable slide clamps, rotating male luer lock, Two Posiflow™ Needleless Valves.	50
--------	--	----

Priming Volume: 4.6mL, Length: 17"



PF-028	Med-Safe™, low volume minibore tubing, quadfurcated set, female luer locks, 0.2 micron filter, non-removable slide clamps, check valves, rotating male luer lock, three Posiflow™ Needleless Valves.	50
--------	--	----




NICU Product Listing

EXTENSION SETS

CATALOG NUMBER	DESCRIPTION	QTY/BOX
----------------	-------------	---------


PF-029	Med-Safe™, low volume minibore tubing, Five-leg set, female luer locks, 0.2 micron filter, non-removable slide clamps, check valves, rotating male luer lock, four Posiflow™ Needleless Valves.	50
--------	---	----

Priming Volume: 4.3mL, Length: 17"




RC-004	Med-Safe™, minibore tubing, quadfurcated set, check valves, female luer lock, non-removable slide clamps, fixed male luer lock.	50
--------	---	----

Priming Volume: 0.8mL, Length: 4"



SM-037	Med-Safe™, minibore tubing, trifurcated set, non-removable slide clamps, rotating male luer lock, three SmartSite® Needleless Valves.	50
--------	---	----

Priming Volume: 0.8mL, Length: 5"




NICU Product Listing

INTOPO™ Enteral Feeding System



INTOPO Syringes: compatible with syringe pumps

CATALOG NUMBER	DESCRIPTION	QTY/BOX
7773	Adapter to convert to oral syringe	100
INT-01	1cc syringe	100
INT-03	3cc syringe	100
INT-05	5cc syringe	100
INT-10	10cc syringe	100
INT-20	20cc syringe	100
INT-30	30cc syringe	50
INT-60	60cc syringe	50



INTOPO Feeding Tubes

CATALOG NUMBER	DESCRIPTION	QTY/BOX
ENPU-20-50	5 French, 20" Polyurethane feeding tube for use up to 30 days	10
ENPU-36-50	5 French, 36" Polyurethane feeding tube for use up to 30 days	10
ENPU-20-65	6.5 French, 20" Polyurethane feeding tube for use up to 30 days	10
ENPU-36-65	6.5 French, 36" Polyurethane feeding tube for use up to 30 days	10
ENPU-20-80	8 French, 20" Polyurethane feeding tube for use up to 30 days	10
ENPU-42-80	8 French, 42" Polyurethane feeding tube for use up to 30 days	10
ENPU-42-100	10 French, 42" Polyurethane feeding tube for use up to 30 days	10
ENPV-16-50	5 French, 16" PVC feeding tube for use up to 3 days	50
ENPV-16-65	6.5 French, 16" PVC feeding tube for use up to 3 days	50
ENPV-16-80	8 French, 16" PVC feeding tube for use up to 3 days	50
ENPV-36-50	5 French, 36" PVC feeding tube for use up to 3 days	50
ENPV-42-80	8 French, 42" PVC feeding tube for use up to 3 days	50

NICU Product Listing

INTOPO™ Enteral Feed-

EXTENSION SETS

CATALOG NUMBER	DESCRIPTION	QTY/ BOX	
EB-024-01	Orange-striped low volume minibore tubing, female luer lock, non-removable slide clamp, universal enteral connector. Priming Volume: 0.50mL, Length: 24"	100	
EB-060-14	Orange-striped low volume minibore tubing, female luer lock, non-removable slide clamp, enteral feeding direction tag, universal enteral connector. Priming Volume: 1.30mL, Length: 60"	100	
EB-060-15	Orange-striped low volume minibore tubing, female luer lock, non-removable slide clamps, feeding tube connector, universal dual enteral connector. Priming Volume: 1.7mL, Length: 60"	100	
EB-060-16	Orange-striped low volume minibore tubing, female luer lock, non-removable slide clamp, male luer lock. Priming Volume: 1.2mL, Length: 60"	100	
EB-060-17	Orange-striped low volume minibore tubing, female luer lock, non-removable slide clamps, feeding tube connector, mini bi-connector w/rotating male luer lock. Priming Volume: 1.4mL, Length: 60"	100	
EB-060-18	Orange-striped low volume mini bore tubing, oral female adapter, non-removable slide clamp, universal enteral connector. Priming Volume: 1.3mL, Length: 60"	100	
EB-060-19	Orange-striped low volume minibore tubing, oral female adapters, non-removable slide clamps, non-vented tethered cap, universal enteral connector. Priming Volume: 1.7mL, Length: 60"	100	
EB-060-20	Orange-striped low volume minibore tubing, oral female adapters, non-removable slide clamps, non-vented tethered cap, mini bi-connector. Priming Volume: 1.4mL, Length: 60"	100	
EB-060-22	Opaque orange minibore tubing, female luer lock, non-removable slide clamp, universal enteral connector. Priming Volume: 2.0mL, Length: 60"	100	
EB-060-23	Opaque orange minibore tubing oral female adapter, non-removable slide clamp, universal enteral connector. Priming Volume: 2.0mL, Length: 60"	100	
EB-060-24	Orange-striped low volume minibore tubing, oral female adapter, non-removable slide clamp, universal enteral connector, INTOPO syringe converter Priming Volume: 1.3 ml, Length: 60"	100	

Acacia™

Ph 714.257.0470 Fx 714.257.0513 • www.acaciainc.com

IMA
Independent Medical Associates

Products distributed by Independent Medical Associates: www.I-MA.com · info@I-MA.com · p: 888.548.4462 · f: 888-548-1462