

IMPROVE
ENTERAL SAFETY

NUTRIO SIMPLEFEED™ INFUSOR

INTRODUCING A DEDICATED ENTERAL FEEDING INFUSOR

Standard electro-mechanical I.V. pumps that have been repurposed for enteral feeding can lead to medication/feeding errors. Acacia Neonatal's NuTrio SimpleFeed Infusor has been designed to look, feel and act differently than electrical pumps in order to eliminate those errors.

It is a safety inspired, neonate-focused, economically sound infusor that offers easier transitions throughout your unit while providing the safest enteral feeding to your patients.

SimpleFeed Infusor tubing is available for the most common delivery rates of 0.5ml/min, 1.0ml/min and 2.0ml/min. They are color coded for easy identification and include our innovative **TwistLok™** technology.

COLOR-CODED TAGS ON THE TUBING

NuTrio SimpleFeed Infusor tubing is calibrated to deliver feeds at 0.5ml, 1.0ml or 2.0ml per minute. Select the correct tubing based on the desired rate of infusion.

Tubing is compatible with Acacia Neonatal 30ml and 60ml syringes.*

DELIVERY RATE:

0.5mL/min
(30mL/hr)

DELIVERY RATE:

1.0mL/min
(60mL/hr)

DELIVERY RATE:

2.0mL/min
(120mL/hr)



*Acacia Neonatal syringes must be used to achieve proper infusion rate.

ORDERING INFORMATION

CAT. NO.	LENGTH / DESCRIPTION / PRIMING VOLUME			
NuTrio SimpleFeed Infusor NTP60				1 BOX
NuTrio SimpleFeed Infusor Winder NTW60				1 BOX
NuTrio SimpleFeed Infusor Tubing Orange-Striped Minibore Tubing, TwistLok Oral Female Adapter, Pinch Clamp, TwistLok Male Adapter				
60ml Syringe Size (Volume To Be Dispensed 20ml - 60ml)				
NTL60-05	60	Delivery Rate 0.5ml/min	0.60ml	50 BOX
NTL60-10	60	Delivery Rate 1.0ml/min	0.95ml	
NTL60-20	60	Delivery Rate 2.0ml/min	1.00ml	
30ml Syringe Size (Volume To Be Dispensed 10ml - 30ml)				
NTL30-05	60	Delivery Rate 0.5ml/min	0.55ml	50 BOX
NTL30-10	60	Delivery Rate 1.0ml/min	0.80ml	
NTL30-20	60	Delivery Rate 2.0ml/min	1.00ml	
NuTrio Syringes				
NTL60	NuTrio 60ml Syringe			50 BOX
NTL30	NuTrio 30ml Syringe			

Contact your IMA representative
or call customer service for more information.

IMA
Independent Medical Associates

11733 66th Street N

Largo, FL 33773

Ph: 888.548.4462

Fax: 888.548.1462

www.I-MA.com

E-Mail: info@I-MA.com

Acacia
Neonatal

www.AcaciaNeonatal.com

© 2012 Acacia, Inc.

Acacia is a registered trademark of Acacia, Inc.

SimpleFeed is a trademark of Acacia, Inc.

M-0104 Rev 040212