

ELIMINATE
MISCONNECTIONS

NUTRIO ENTERAL FEEDING SYSTEM™

GRAVIFEED™: THE WORLD'S FIRST GRAVITY FEEDING SPECIFIC SYRINGE

Eliminate waste and save money while improving safety with the GraviFeed gravity feeding syringe.

- Built in TwistLok™ feature prevents accidental disengagement and dangerous misconnections between the feeding tube and GraviFeed syringe.
- Snap tight sealable lid prevents exposure to airborne particulates or contamination and the unnecessary loss of precious milk or formula.
- Reusable primer effortlessly creates positive pressure and eliminates the need for a plunger, which saves the hospital money.
- An easier way to deliver gravity feeding without the push and pull associated with a traditional plunger.
- Hydrophobic filter built into the lid allows air to pass through without letting enteral fluid leak out, making the GraviFeed a perfect solution for post feed venting.
- Non-DEHP, Latex Free.



NUTRIO ENTERAL FEEDING SYSTEM



Traditional
Venting




GraviFeed
Venting



Perfect for Venting – Venting with GraviFeed is as simple as snapping the lid closed. No gauze or hanging necessary.



Less Waste – GraviFeed utilizes a reusable primer to create positive pressure to initiate flow. This method eliminates waste because an expensive single use syringe plunger is not being thrown away with every feed. An optional single use plunger is available for hospitals not interested in switching to the reusable primer.

ORDERING INFORMATION

| CAT. NO. | DESCRIPTION | |
|---|---|-----------|
|  | | 50 BOX |
| GF-30 GraviFeed | 30ml Syringe | |
|  | | 50 BOX |
| GF-60 GraviFeed | 60ml Syringe | |
|  | | 10 BOX |
| GF-001 GraviFeed | Plunger (For use with GF-30 and GF-60) | |

| CAT. NO. | DESCRIPTION | |
|--|------------------------|-----------|
|  | | 10 BOX |
| GF-002 GraviFeed | Primer | |
| GF-002L | Large GraviFeed Primer | |
|  | | 1 EACH |
| GF-003 GraviFeed | Rack | |

Distributed By:



11733 66th St N Suite 113
Largo, FL 33773

PH: 888.548.4462
FX: 888.548.1462

E-Mail: info@I-MA.com
Web: www.I-MA.com

© 2011 Acacia, Inc.
Acacia is a registered trademark of Acacia, Inc.
TwistLok, GraviFeed and NuTrio Enteral Feeding System are trademarks of Acacia, Inc.
Patent Pending