

# Posiflow™ Valve Product Listing

The Posiflow™ needleless valve is engineered to automatically reverse blood reflux associated with conventional needleless valves. This positive fluid displacement feature will ensure the patency of your catheters which is particularly important in the neonatal population. The luer activated straight through fluid path reduces the chance of hemolysis and allows for easy administration of fluids. The Posiflow™ valve has a flat swabbable surface that reduces the chance of microbial agents entering the fluid path. It is compatible for use with lipids and chemotherapeutic agents.

## FEATURES:


- Latex-Free
- Length: 1.25"
- Priming Volume: 0.06 ml
- Residual Volume: 0.03 ml
- Displaced Volume: 0.03 ml
- Maximum Flow Rate<sup>1</sup>: 3.2 L/hr



CATALOG NUMBER	DESCRIPTION	QTY/BOX
----------------	-------------	---------

PF-001	Minibore tubing, non-removable slide clamp, male luer lock, Posiflow™ needleless valve	50
--------	--	----

Priming Volume: 0.4mL, Length: 10"




PF-002	Standard bore tubing, non-removable slide clamp, male luer lock, removable Posiflow™ needleless valve	50
--------	---	----

Priming Volume: 1.0mL, Length: 10"




PF-003	Minibore tubing, non-removable slide clamp, male luer lock, Posiflow™ needleless valve	50
--------	--	----

Priming Volume: 0.53 mL, Length: 14"



PF-004	Minibore tubing, non-removable slide clamp, Posiflow™ needleless valve t-connector, female luer lock	50
--------	--	----








Priming Volume: 1.5 mL, Length: 38"










Acacia™

Ph 714.257.0470 Fx 714.257.0513 • www.acaciainc.com

# Posiflow™ Valve Product Listing

CATALOG NUMBER	DESCRIPTION	QTY/BOX	
PF-005	<p>Macrobore tubing, 20 drops/ml universal drip chamber, roller clamp, male luer lock, Posiflow™ needleless valve</p> <p>Length: 100"</p>	50	
PF-006	<p>Macrobore tubing, 20 drops/ml universal drip chamber, roller 0.2 micron filter, male luer lock, Posiflow™ needleless valve</p> <p>Length: 100"</p>	50	
PF-007	<p>22 mm vial adapter with Posiflow™ needleless valve</p>	50	
PF-008	<p>Trifurcated set, low volume minibore tubing, female luer locks, 0.2 micron filter, non-removable slide clamps, male luer lock, Posiflow™ needleless valve</p> <p>Priming Volume: 3.6 mL, Length: 15"</p>	50	
PF-009	<p>Trifurcated set, low volume minibore tubing, female luer locks, non-removable slide clamps, male luer lock, Posiflow™ needleless valve</p> <p>Priming Volume: 0.8 mL, Length: 14"</p>	50	
PF-010	<p>Minibore tubing, female luer lock, non-removable slide clamps, mini bi-connector with luer lock, Posiflow™ needleless valve</p> <p>Priming Volume: 0.7ml, Length: 11"</p>	50	
PF-011	<p>Macrobore tubing, 20 drops/ml universal drip chamber, back check valve, roller clamp, 1.2 micron filter, rotating male luer lock, removable back check valve, Posiflow™ needleless valve</p> <p>Length: 100"</p>	50	

# Posiflow™ Valve Product Listing

CATALOG NUMBER	DESCRIPTION	QTY/BOX	
PF-013	Minibore tubing, female luer lock, non-removable slide clamp, male luer lock, removable Posiflow™ needleless valve  Priming Volume: 0.4 mL, Length: 11"	50	
PF-014	22 mm vial adapter	50	
PF-015	Minibore tubing, non-removable slide clamp, rotating male luer lock, Posiflow™ needleless valve  Priming Volume: 0.4 mL, Length: 10"	50	
PF-016	Low volume macrobore tubing, female luer lock, pinch clamp, rotating male luer lock, removable Posiflow™ needleless valve  Priming Volume: 1.0 mL, Length: 8"	50	
PF-017	Minibore tubing, t-injection site with rotating male luer lock, Posiflow™ needleless valve (BTO)  Priming Volume: 0.25 mL, Length: 5.5"	50	
PF-018	Minibore tubing, t-injection site with male luer slip, Posiflow™ needleless valve  Priming Volume: 0.3 mL, Length: 6"	50	
PF-019	Trifurcated set, low volume minibore tubing, female luer locks, 0.2 micron filter, non-removable slide clamps, rotating male luer lock, Posiflow™ needleless valve  Priming Volume: 3.6 mL, Length: 16"	50	


Acacia™

Ph 714.257.0470 Fx 714.257.0513 • [www.acaciainc.com](http://www.acaciainc.com)

# Posiflow™ Valve Product Listing

**CATALOG NUMBER**      **DESCRIPTION**      **QTY/BOX**

PF-020      Macrobore tubing, female luer lock, large non-removable slide clamp, Posiflow™ needleless valve y-connector, rotating male luer lock




Priming Volume: 1.25 mL, Length: 6"

PF-021      Minibore tubing, non-removable slide clamp, rotating male luer lock, Posiflow™ needleless valve




Priming Volume: 0.15 mL, Length: 4"

PF-022      Quadfurcated set, minibore tubing, female luer locks, non-removable slide clamps, rotating male luer lock, Posiflow™ needleless valve




Priming Volume: 0.6 mL, Length: 4"

PF-023      Trifurcated set, low volume minibore tubing, female luer locks, non-removable slide clamps, rotating male luer lock, Posiflow™ needleless valve




Priming Volume: 0.5 mL, Length: 4"

PF-024      Trifurcated set, minibore tubing, female luer lock, 0.2 micron filter, back check valves, non-removable slide clamp, rotating male luer lock, Posiflow™ needleless valve



Priming Volume: 4.6 mL, Length: 17"

PF-025      Bifurcated set, minibore tubing, female luer lock, non-removable slide clamps, rotating male luer lock, Posiflow™ needleless valve









Priming Volume: 0.5 mL, Length: 5"

Acacia™

Ph 714.257.0470 Fx 714.257.0513 • www.acaciainc.com

# Posiflow™ Valve Product Listing

CATALOG NUMBER	DESCRIPTION	QTY/BOX	
PF-026	Trifurcated set, low volume minibore tubing, female luer locks, non-removable slide clamps, rotating male luer lock, Posiflow™ needleless valves	50	 <p>Priming Volume: 0.7 mL, Length: 5"</p>
PF-027	Trifurcated set, minibore tubing, female luer lock, 0.2 micron filter, back check valve, non-removable slide clamps, rotating male luer lock, Posiflow™ needleless valves	50	 <p>Priming Volume: 4.6 mL, Length: 17"</p>
PF-028	Quadfurcated set, low volume minibore tubing, female luer locks, 0.2 micron filter, non-removable slide clamps, back check valves, rotating male luer lock, Posiflow™ needleless valves	50	 <p>Priming Volume: 4.2 mL, Length: 17"</p>
PF-029	Quintfurcated set, low volume minibore tubing, female luer locks, 0.2 micron filter, non-removable slide clamps, back check valves, rotating male luer lock, Posiflow™ needleless valves	50	 <p>Priming Volume: 4.3 mL, Length: 17"</p>
PF-030	Minibore tubing, non-removable slide clamp, rotating male luer lock, Posiflow™ needleless valve	50	 <p>Priming Volume: 0.2 mL, Length: 7"</p>
PF-031	Minibore tubing, non-removable slide clamp, rotating male luer lock, removable Posiflow™ needleless valve	50	 <p>Priming Volume: 0.35 mL, Length: 9"</p>


Acacia™

Ph 714.257.0470 Fx 714.257.0513 • www.acaciainc.com


# Posiflow™ Valve Product Listing

CATALOG NUMBER	DESCRIPTION	QTY/BOX
----------------	-------------	---------

PF-032	Minibore tubing, female luer lock, 0.2 micron filter, non-removable slide clamps, rotating male luer lock, Posiflow™ needleless valve  Priming Volume: 4.2 mL, Length: 10"	50
--------	--	----



PF-033	Trifurcated set, low volume minibore tubing, female luer locks, non-removable slide clamps, rotating male luer lock, Posiflow™ needleless valve  Priming Volume: 0.9 mL, Length: 8"	50
--------	---	----



Acacia™

Ph 714.257.0470 Fx 714.257.0513 • [www.acaciainc.com](http://www.acaciainc.com)

**IMA**  
Independent Medical Associates

Products distributed by Independent Medical Associates: [www.I-MA.com](http://www.I-MA.com) · [info@I-MA.com](mailto:info@I-MA.com) · p: 888.548.4462 · f: 888-548-1462