

**REDUCE  
INFECTION**

# **SAFESAMPLE™ BLOOD SAMPLING SET**



## **SAFESAMPLE™ UAC BLOOD GAS SAMPLING**

In the critically ill newborn, frequent blood sampling from an umbilical artery catheter (UAC) may be necessary. These fragile newborns, whether born pre-mature or ill, require the gentlest of care particularly in preserving their blood. The SafeSample™ Blood System allows for frequent blood gas sampling from a UAC, without unnecessary trauma to the newborn.

In comparison to using a traditional open stopcock for sampling, the system's closed valves significantly reduce the "residual volume" in the valve and further contribute to preventing bacterial colonization. The system also reduces blood exposure to the clinician by using closed needleless valves.

# SAFESAMPLE BLOOD SAMPLING SET

## PRODUCT FEATURES & BENEFITS:

- Sampling port and clearance port both use closed needleless valves
- Automatic clearing of valve when flushing eliminates the need to manipulate the valve handle
- No cannulas needed
- Reduces blood exposure to clinicians
- Closed design reduces the risk of microbial colonization
- Color-coded valve indicates sampling port
- Universal system can be used with most transducers
- Easy to prime and use

ORDERING INFORMATION

### SS-001

Blood sampling system for use with prefilled flush syringes. Accomodates clearance draws of 1.0-2.0ml.

TOTAL APPROX. PRIMING VOLUME: 3.8ml

APPROX. LENGTH: 9

25  
BOX



**IMA**  
Independent Medical Associates

11733 66th Street North  
Suite 113, Largo, FL 33773

PH: 888.548.4462

FX: 888.548.1462

E-Mail: [info@I-MA.com](mailto:info@I-MA.com)

Web: [www.I-MA.com](http://www.I-MA.com)

© 2011 Acacia, Inc.  
SafeSample is a registered trademark of Acacia, Inc.

M-0112  
071111