


# MedSafeDemo

## Standard Multi-Access Extension Set

Quadfurcated extension set, minibore tubing, female luer locks, non-removable slide clamps, rotating male luer lock


Approx. Priming Volume: Varies    Approx. Length: 7"



## Acacia's MedSafe Extension Set #BG-005-01

Quadfurcated extension set, minibore tubing, check valves, female luer lock, non-removable slide clamps, rotating male luer lock

Approx. Priming Volume: 0.65ml    Approx. Length: 5"



\* Use Fine Point Sharpie

Acacia®

Ph 714.257.0470

Fx 714.257.0513

www.acaciainc.com

**IMA**  
Independent Medical Associates

Products distributed by Independent Medical Associates: [www.I-MA.com](http://www.I-MA.com) · [info@I-MA.com](mailto:info@I-MA.com) · p: 888.548.4462 · f: 888-548-1462