

# Subcutaneous Needle Infusion Set

Hypodermoclysis, the subcutaneous infusion of fluids, is a useful and easy hydration technique suitable for mildly to moderately dehydrated adult patients, especially the elderly. The method is considered safe and does not pose any serious complications. The most frequent adverse effect is mild subcutaneous edema that can be treated by local massage or systemic diuretics. Approximately 3 L can be given in a 24-hour period at two separate sites. Common infusion sites are the chest, abdomen, thighs and upper arms.



Designed specifically for the delivery of medications directly to the subcutaneous tissue.

Catalog number: MH4000  
Box quantity: 50  
Approximate length: 42"  
Approximate priming volume: 0.25 ml

Microbore tubing with 27g x 1/2" white wing needle

Catalog number: MH4001  
Box quantity: 50  
Approximate length: 24"  
Approximate priming volume: 0.17 ml

Microbore tubing with 27g x 1/2" white wing needle

Acacia™

Ph 714.257.0470 Fx 714.257.0513 • [www.acaciainc.com](http://www.acaciainc.com)

**IMA**  
Independent Medical Associates

Products distributed by Independent Medical Associates: [www.I-MA.com](http://www.I-MA.com) · [info@I-MA.com](mailto:info@I-MA.com) · p: 888.548.4462 · f: 888-548-1462