



# THE BLOCKS™

Radio Translucent Pre-Wired Electrodes

## Radiolucent:

Radiolucent electrode wires provide better visulation for x-ray

## Improved Flexibility:

Softer lead wire improves flexibility with additional length

## Better Fit:

Molded tapered connector for better fit into patient cable

## Better Signal:

Center lead wire allows additional strenght and provides a better ECG signal

## Superior Hydrogel:

Provides excellent attachment with no residue when removed or repositioned

## Available in 3 Sizes:

Premie (7/8 x 7/8) - Order # PC1-6C - (1 Bx 100ct)

Newborn (7/8 x 1-1/4) - Order # PC2-6C - (1 Bx 100ct)

Pediatric (1-1/4 x 1-1/4) - Order # PC3-6C - (1 Bx 100ct)

Distributed by:

**IMA**  
Independent Medical Associates

Ph. 888-548-4462 Fax. 888-548-1462