

CharterMedical


A Lydall Company



Just Say NO to Needles!

The Charter Medical NAIL (Needle-less Access Injection Line) is a unique alternative to existing sampling site couplers and potentially hazardous needle/syringe systems. Just spike the set into any solution or blood bag and attach a syringe to the female luer lock (syringe adapter). It's safe, efficient and cost-effective.

- Helps comply with blood-borne pathogen safe-needle laws nationwide
- Permits attachment to a flexible non-vented container to prepare or process blood components
- For use in preparation of RBC, platelet and cryoprecipitate aliquots
- Sterile connect for quality control sample collection
- Many applications for use in blood component transfer:
 - Use as a needle-less sampling port
 - Access autologous blood for platelet count or lab testing
- Customer-requested product
- Latex-free



CATALOG #	DESCRIPTION	CASE QTY
03-220-50	3", pinch clamp with piercing pin, dual needle adapter and spike port	96
03-220-90	4", pinch clamp with piercing pin and female luer lock	96
03-220-92	2", pinch clamp with piercing pin and female luer lock	96
03-220-94	2", pinch clamp with piercing pin and male luer lock	96
03-220-95	2", pinch clamp with male luer lock and female luer lock	96

Distributed by:

Independent Medical Associates

IMA

Independent Medical Associates

11733 66th Street North
Suite 113
Largo, FL 33773

Phone: 888-548-4462

Fax: 888-548-1462

E-Mail: Info@I-MA.com

The Federal Needlestick Safety and Prevention Act (H.R. 5178) mandates needlestick safety practices nationwide. A copy of the document is available from Charter Medical Ltd. Customer Service.

Charter Medical, Ltd.
A Lydall Company
3948-A Westpoint Blvd.
Winston-Salem, NC 27103
Phone: 336.768.6447
Fax: 336.714.4241
Toll Free: 1.866.458.3116

08-05 NAIL (7,500)