

Platelet Sampling Devices

For use with bioMérieux BacT/ALERT® Bacterial Detection Systems

Our easy-to-use dock, fill and inoculate method transfers platelet samples safely and efficiently. These customer-developed, widely-utilized systems feature:

- More syringe capacity than previously available sampling devices
- Integrated bioMérieux culture vial adapter cap
- Syringe configurations for accurate sample measurement
- Tubing compatible with sterile connection devices



Catalog #	Description	Case QTY
03-220-TM	Transfer Set with 100 mL Bag and Safety Needle Assembly	48
03-220-BC	Transfer Set with 10 mL BD™ Syringe and Cap/Vial Adapter	24
03-220-BC12	Transfer Set with 12 mL Syringe and Cap/Vial Adapter	24
03-220-BC20	Transfer Set with 20 mL BD™ Syringe and Cap/Vial Adapter	24
03-220-TMBC	Transfer Set with 100 mL Bag, 12 mL Syringe and Cap/Vial Adapter	24

BD is a trademark of Becton Dickinson and Co.
References and technical data available upon request.

Distributed by:

CharterMedical

Your Neonatal Transfusion Expert

Charter Medical, Ltd.
A Lydall Company
3948-A Westpoint Blvd.
Winston-Salem, NC 27103
Customer Service 1.866.458.3116
www.chartermedical.com

Independent Medical Associates

IMA
Independent Medical Associates

11733 66th Street North
Suite 113
Largo, FL 33773

Phone: 888-548-4462
Fax: 888-548-1462
E-Mail: Info@I-MA.com