



## BabyPod II™

Transporting nature's precious cargo

### Features

- Cushioned and insulated to protect baby during transport
- Lightweight unit for one-operator transport
- Adaptable to many environments



# BabyPod™ II is designed to safely convey a pediatric patient between departments within the hospital or between medical facilities.

## Protects Baby

- Provides thermo-stable environment.
- Shock absorbent foam absorbs vibration during transport.
- Vacuum mattress conforms to baby's shape to limit unwanted movement.
- Viewing lid protects baby from loose moving objects while allowing visual monitoring.

## Easy to Use

- Lightweight carbon fiber shell
- 20-pound unit suitable for one-operator transport
- Straps quickly to any stretcher or gurney

## Flexible

- Non-electric heat source: TransWarmer® Infant Transport Mattress.
- Infusion pumps and/or IV's attach directly to the unit.
- BabyPod™ II may be hand-carried over rough terrain or when a stretcher is unavailable.

## Multiple uses

- Evacuation
- Search and rescue
- Ground ambulance
- Between hospital transport
- Inter-department transport



**TransWarmer® mattress keeps infants warm while in the BabyPod™ II or, as a stand alone heat source in the NICU, nursery, or anywhere newborns can experience cold stress.**



*Only nature builds a better nest*



## **TransWarmer® Infant Transport Mattress**

TransWarmer® Infant Transport Mattress provides a thermo-stable environment during transport for up to two hours. The TransWarmer Infant Transport Mattress is a gel-filled, disposable, thermal mattress designed for use in neonatal transports. TransWarmer's patented WarmGel® filler provides even heating to minimize the risk of hypothermia.

### **Features & Benefits**

*Minimizes the risk of cold stress*

- Patented WarmGel® mattress
  - Nests infant, absorbing vibration during transport
- Consistent, uniform heat distribution
  - Prevents the complications of cold stress
- Heats in less than 60 seconds
  - Provides rapid thermal support
- Easy to activate
  - Simple, no inner bags to rupture
- Non-toxic, food-grade gel

# BabyPod II™

Transporting nature's precious cargo

## Infant Transport Device

BabyPod™ II consists of a lightweight carbon fiber outer shell, lined with a shock absorbent foam inner layer, and has transparent shields for viewing the baby. It also contains a patient positioning vacuum mattress, stretcher fixing straps to secure the device for transport, and safety straps to secure the infant inside the device during transport.



## Vital Statistics

Device Dimensions	38" x 17" x 12"
Device Weight (No accessories)	20lbs/9kg
Patient Weight Maximum	Up to 20lbs/9kg

## Ordering Information

70002	BabyPod™ II Infant Transport Device (1 each)
20421	TransWarmer® Infant Transport Mattress (6 per box)

BabyPod II™



**IMA**  
Independent Medical Associates

Keeping You at the Forefront  
of Women's Health Care™

Transwarmer® is a registered trademark of Prism Enterprises, L.P.

11733 66th St. N #113 Largo, FL 33771 • 888.548.4462 • 727.548.4462

www.I-MA.com •

©2006 CooperSurgical, Inc.

Transwarmer and WarmGel are registered trademarks of Prism Enterprises, a CooperSurgical company.

Form # 80641 Rev. 08/06