



Belly Bands™

Abdominal Fetal Monitoring Belt

Features

- Flexible and easy to use — allows mother to place band comfortably over abdomen
- Comes in precut versions or in a 10 yard roll
- 14 inch pre-cut length — covers the belly and secures fetal monitor placement
- Breathable material - helps wick away moisture



Belly Bands™

Abdominal Fetal Monitoring Belt



Instructions

- Choose appropriate size
- Open Belly Band and gather with both hands
- Slip Belly Band over both feet
- Slip Belly Band comfortably over abdomen



Ordering Information

13010	10" Width Belly Bands on Roll (10"x10yds.)	1 Roll
13011	10" Width Belly Bands - Pre-cut (10"x14")	25 per Box
13012	12" Width Bell Bands on Roll (12"x10yds.)	1 Roll
13013	12" Width Belly Bands - Pre-cut (12"x14")	25 per Box
13014	14" Width Belly Bands on Roll (14"x10yds.)	1 Roll
13015	14" Width Belly Bands - Pre-cut (14"x14")	25 per Box
13016	16" Width Belly Bands on Roll (16"x10yds.)	1 Roll
13017	16" Width Belly Bands - Pre-cut (16"x14")	25 per Box
13019	18" Width Belly Bands - Pre-cut (18"x14")	25 per Box

IMA
Independent Medical Associates

11733 66th St. N #113 Largo, FL 33771 • 888.548.4462 • 727.548.4462
Fax 888.548.1462 • www.I-MA.com

Form # 81238 Rev. 06/07
©2007 CooperSurgical, Inc.

Belly Bands™ is a trademark of CooperSurgical.