



**EZ-HSG™**

**Simple, cost effective choice for  
ultrasound examination of the uterus**

**Features**

- Balloon-less
- Creating a single channel for easier maneuverability
- "Soft Foam" stopper replaces balloon
- Simplifies saline-infused procedures
- Latex-free





Simple, cost effective choice for ultrasound examination of the uterus

### 3 Different Sizes

- Small – for nulligravida and post-menopausal patients
- Medium – for primipara and multipara patients
- Large – for incompetent cervix

EZ-HSG is designed for Sonohysterography along with Hysterosalpingography. Simply insert the catheter's soft foam stopper into the internal canal. Slowly inject fluid to distend the uterine cavity.



### Dimensions

Small SM4200	Catheter O.D. 1.6 mm External cannula 3 mm Foam cone 6 mm flared to 12 mm
Medium SM4400	Catheter O.D. 2 mm External cannula 3 mm Foam cone 10 mm flared to 15 mm
Large SM4600	Catheter O.D. 3 mm External cannula 5 mm Foam cone 8 mm flared to 18 mm
All catheters	Overall length 35 mm Catheter tip length 1 to 1.5 cm Foam length 2.5 to 3 cm

Distributed By:



Independent Medical Associates

11733 66th Street North  
Suite 113  
Largo, FL 33773

Phone: 888-548-4462  
Fax: 888-548-1462  
E-Mail: Info@I-MA.com