

Hysterosonography



GYNECATH

Uterine Infusion Catheter

Comfort and control in a single device

Features

- Comfortable
 - Soft, flexible catheter minimizes discomfort during uterine exams
 - Gentle on delicate tissues, flexible plastic catheter and latex-free balloon
 - Minimally invasive
- Easy to control
 - Physician can control catheter with one hand
 - Once catheter is placed, there is no need for further manipulation
 - Significantly reduces the need for metal tenaculum
- Easy to use
 - Sliding plastic sheath aids in catheter placement
 - Latex-free balloon inflates easily once catheter is inserted
 - Balloon reduces fluid leakage, allowing optimal diagnostic viewing
 - Available in a procedure tray set to minimize exam time



GYNECATH

Comfort and control in a single device

GyneCath

The GyneCath Uterine Infusion Catheter is a minimally invasive, easy-to-use device that ensures a diagnostically enhanced, yet gentle, patient exam. Designed for use in saline-infusion hysterosonography and water-based hysterosalpingography.

Enhanced Diagnostic Images

GyneCath is designed for use in saline-infusion transvaginal ultrasound and in radiologic exams using water-based contrast media. Once the balloon is inflated and fluid infused, the clinician's view of the uterine cavity and fallopian tubes is significantly enhanced.

Product Specifications

Radiopaque 12" flexible 5 French/ 7 French catheter has a rounded closed distal tip with a side port proximal to the tip. The attached latex-free balloon has a 3cc/6cc capacity and is proximal to the side port. The dual lumen tube allows for balloon inflation and injection of the contrast media or saline solution through the side port.

GyneCath Uterine Infusion Catheter

- Two sizes available:
 - 5 French for non-parous patients
 - 7 French for previously pregnant patients
- Available in tray sets for clinician convenience.



Procedure Tray Sets

Catheter, fenestrated drape, disposable speculum, swabs and swab sticks, sterile wraps, 4" x 4" gauze pad, sanitary napkin, lubricant, providone-iodine solution, 36" extension tube, and 12cc syringe/18 gauge needle packed in a sterile tray.

Order Information

Catheters (10 per box 40 per case)

GyneCath
19610 5F catheter
19611 7F catheter

Tray Sets (10 trays per case)

GyneCath
19612 5F catheter
19613 7F catheter
Ob/Gyn
19614 Tray set (no catheter)

Distributed By:

IMA
Independent Medical Associates

Independent Medical Associates

11733 66th Street North
Suite 113
Largo, FL 33773

Phone: 888-548-4462
Fax: 888-548-1462
E-Mail: Info@I-MA.com