


# INTRODUCING Mystic<sup>®</sup> II



Enhancing a winner —  
The Mityvac line  
of vacuum assisted  
delivery systems

*Remember, one size cup does  
not fit all deliveries!*

**IMA**  
Independent Medical Associates

11733 66th St. N #113 Largo, FL 33771 • 888.548.4462 • www.I-MA.com  
©2007 CooperSurgical, Inc.

## Safety

Lower scalp lacerations by 68%

## Low Profile

Mushroom<sup>®</sup> Cup flexes up to 90° in any  
direction for easier insertion

## Superior Control

Handle can rotate 180° without transferring  
rotation to the cup

## Versatile

An array of cups to choose from

Mystic<sup>®</sup> is a product of CooperSurgical, Inc.

Form # 81237OBYN 06/07