



NEO-fit™

For easy, secure endotracheal tube stabilization

Key Features & Benefits

- Provides a sure means of securing the tube to the patient
- Unique design eliminates the use of tape on the endotracheal tube and minimizes the use of tape on the patient
- Easy to use; "one size fits all"
- Demarcations are clearly seen on tube
- Hook and loop strap allows rapid release and re-fit making tube adjustment simple
- Longer term use and cost efficient



NEO-fit™ - For easy, secure endotracheal tube stabilization

The NEO-fit provides quick, secure endotracheal tube stabilization. Designed to firmly grip a Neonatal endotracheal tube, NEO-fit utilizes metal clips to mechanically grip the soft material of the tube. The device is positioned on baby's face using a hydrocolloid adhesive that is hypoallergenic and velcro strap allows for rapid release and readjustment of the tube.

Ordering information

NEO-fit includes one plastic tube grip with adjustable hook and loop strap, two adhesive foam pads and two foam "lollypop" securing tapes (optional use).

Product Number
42-2540

Packaging

Non-sterile, 20per box, 80 per case (4 boxes)

Complimentary Products:

NEO-prep - Provides protection for the tender neonatal septum

NEO-prep provides extra skin barrier, minimizing septal breakdown and can be used with CooperSurgical's NEO-Fit and INCA products. The pre-cut hypoallergenic skin barrier is available in either Small or Large size.

Ordering information

Product Number	Size
44-8000	Small
44-8001	Large

Packaging

8 per strip, 10 strips per box, 3 boxes (240) per case



www.I-M-A.com

11733 66th St. N Suite 113 Largo FL 33771
888.548.4462 • 727.548.4462