


Stainless Steel Instruments for Gynecologic Procedures




Vu-More Wide Rear-Opening Speculum

- Permits easy removal from the patient without disturbing the catheter or instrument
- The longer blade makes it suitable for obese patients

Product No. 63-902 Medium

Product No. 63-900 Large



Weisman-Graves Open-Sided Speculum


- Speculum design permits easy removal without disturbing catheter or instruments
- Ideal time-saving device for Hysterosonography
- Easily adjustable to most anatomic needs
- Available with either left or right side opening

Product No. 391-667A Left side open – Medium

Product No. 391-668A Right side open – Medium

Product No. 391-667 Left side open – Large

Product No. 391-668 Right side open – Large



Standard Luer Needle Extension

- Allows use of less expensive standard size needles for paracervical blocks
- Malleable material can be bent to angle needed



Product No. 433-131A Reusable 3" extension

Product No. 433-132 Reusable 4" extension


Packaging

All of the products listed above are packaged individually in a clear plastic bag. Four per box. Non-sterile.

Goldstein GRASP™ Cervical Stabilizer



Normal tenaculum is designed for deep penetration which often causes considerable pain and bleeding.



GRASP provides thin pincers which are designed to stabilize the uterus with minimal pain and bleeding.

The **GRASP** provides rapid and relatively atraumatic stabilization of the cervical os during placement of various catheters, endometrial samplers and other intrauterine devices.

- Thin pincers provide stabilization with minimal pain and bleeding.
- Single-click ratchet prevents over-tightening of device.
- Reusable stainless steel can be steam autoclaved.

Ordering Information

Goldstein GRASP Cervical Stabilizer

Product No. 64-6026

Packaging — One each, non-sterile.

Distributed By:

IMA
Independent Medical Associates

Independent Medical Associates

11733 66th Street North
Suite 113
Largo, FL 33773

Phone: 888-548-4462
Fax: 888-548-1462
E-Mail: Info@I-MA.com