



TAMPA CATHETER

The simple, cost-effective choice for
ultrasound examination of the uterus

Features

- Unique introducer sheath facilitates insertion into cervical os with minimal trauma.
- 7cm marker provides optimal placement of catheter tip in uterus.
- Balloonless configuration is ideal for post-menopausal patients with stenotic cervical canals, as well as for many nulliparous patients.
- Small catheter size allows easy placement next to ultrasound probe.




TAMPA CATHETER

The simple, cost-effective choice for
ultrasound examination of the uterus

Ordering information

Tampa Catheter

- **PRODUCT NO.** 61-2005
5 French, clear, flexible catheter with
introducer sheath, 33cm long.
- **PACKAGING**
Sterile, 20 sets per box
4 boxes (80 sets) per case



Stainless Steel Instruments for Gynecologic Procedures

Vu-More™ Wide Rear-Opening Speculum


- Permits easy removal from the patient without disturbing the catheter or instrument
- The longer blade makes it suitable for obese patients

Product No. 63-902

Medium

Product No. 63-900

Large



Standard Luer Needle Extension

- Allows use of less expensive standard size needles for paracervical blocks
- Material can be bent to angle needed

Product No. 433-131A

Reusable 3" extension

Product No. 433-132

Reusable 4" extension



Weisman-Graves Open-Sided Speculum

- Speculum design permits easy removal without disturbing catheter or instruments
- Ideal time-saving device for Hysterosonography
- Easily adjustable to most anatomic needs
- Available with either left or right side opening

Product No. 391-667A

Left side open – Medium

Product No. 391-668A

Right side open – Medium

Product No. 391-667

Left side open – Large

Product No. 391-668

Right side open – Large



Packaging

All of the products listed above are packaged individually in a clear plastic bag. Four per box. Non-sterile.

TAMPA CATHETER

Distributed By:

IMA
Independent Medical Associates

Independent Medical Associates

11733 66th Street North
Suite 113
Largo, FL 33773

Phone: 888-548-4462
Fax: 888-548-1462
E-Mail: Info@I-MA.com