




TransWarmer

Only nature builds a better nest

Features & Benefits

Minimizes the risk of cold stress

- Patented WarmGel® mattress
 - Nests infant, absorbing vibration during transport
- Consistent, uniform heat distribution
 - Prevents the complications of cold stress
- Heats in less than 60 seconds
 - Provides rapid thermal support
- Easy to activate
 - Simple, no inner bag to rupture



Transwarmer: It's not just for transport anymore

TransWarmer

Only nature builds a better nest

Product Overview

Transwarmer Infant Transport Mattress is a gel-filled, disposable thermal mattress designed for use in all cases when cold stress is a concern. Patented WarmGel heats quickly to minimize the risk of hypothermia.

Applications

- Prevent or Alleviate Cold Stress
- Neonatal Intensive Care
- Labor and Delivery Suite
- Labor and Delivery Operating Room
- Pediatrics
- Transport
- EMS

Ordering Information

Reference 20421

TransWarmer Infant Transport Mattress
(with patented WarmGel filler)

6 per case



Safe & Convenient

- Non-toxic food-grade gel
- Single-patient use minimizes risk of infection
- Portable heat source requires no electricity

Reliable

- Reaches maximum temperature in less than 60 seconds
- Trigger disc ensures activation when needed

Clinically Proven

- Reduces hypothermia risk in neonates
- Maintains 104° F for up to two hours

The Unique Benefits of WarmGel

- Better heat distribution increases blood flow
- Optimizes temperature maintenance in neonates and small pediatric patients
- Absorbs vibration during infant transport
- Patented gel and soft non-woven surface gently conform to patient's shape for better nesting
- Reduces the risk of cold stress



IMA
Independent Medical Associates

Keeping You at the Forefront
of Women's Health Care

11733 66th St. N #113 Largo, FL 33771 • 888.548.4462 • 727.548.4462

www.I-MA.com •

©2005 CooperSurgical, Inc.

Transwarmer and WarmGel are registered trademarks of Prism Enterprises, a CooperSurgical company.

Form # 80836 Rev. A 01/05