



Diros Technology Inc.




RF Pain Management
& Neurosurgical Solutions




MULTI-LESION ADAPTOR — MLA-4

dirostech.com


Multi Channel Stimulation





Multi Lesion Standard RF



Multi Lesion Pulsed RF



Probe Test



- ✓ Compatible with OWL URF-3AP RF Generator
- ✓ User-Friendly true plug & play system
- ✓ Simplified Probe ID and Probe Management System
- ✓ Ability to create up to 4 Lesions Simultaneously
- ✓ Real-Time Graphical Display for all 4 Lesions
- ✓ Continuous Impedance Monitoring for all 4 Channels
- ✓ Auto/Manual Extended Ramp available for all Lesions
- ✓ Numerous built in visible/invisible Safety Parameters
- ✓ Unique Probe Test Adaptor for all 4 Channels

Manufactured By: DIROS TECHNOLOGY INC.

Distributed By: