




Diros Technology Inc.

OWL[®] RF Pain Management & Neurosurgical Solutions



Diros OWL[®] Sterile Single-Use RF Cannulae

- Lengths Available: 50mm, 54mm, 100mm, 145mm, 175mm
- Gauges Available: 14Ga to 22Ga
- Bare Tip Lengths: 2mm to 15mm
- Bare Tip Configurations: Sharp, Blunt, Straight, Curved



Diros OWL[®] Sterile Single-Use & Re-Usable RF Probes


- Lengths Available: 5cm, 10cm, 15cm, 20cm
- Colour coded for easy reference
- Thermocouple & Thermistor Technology
- Compatible with competitive RF Generators

Diros OWL[®] GD-pad




Diros OWL[®] SMART Cables





Diros RF Pain Management Sets

- Facet Denervation Sets
(Available in 5cm, 10cm, 15cm and 20cm Lengths)
- Facet Rhizotomy Set
- Custom Sets available upon request



Diros RF Neurosurgery Sets

- OWL Tasker Intracranial Set
- OWL Gasserian Ganglion Set
- OWL D.R.E.Z. Set
- OWL Cordotomy Set
- Custom Sets available upon request

Diros URF-3AP RF Generator



MANUFACTURED BY: **DIROS TECHNOLOGY INC.**

DISTRIBUTED BY:



232 Hood Road
Markham, Ontario
L3R 3K8 Canada
Tel: (905) 415-3440
Fax: (905) 415-0667
E-Mail: sales@dirostech.com



Registered Trademark of Diros Technology Inc.

Diros Technology Inc. reserves the right to change specifications or to incorporate design changes without further notice and without incurring any obligations relating to equipment previously manufactured or delivered.

© 2009 Diros Technology Inc. All Rights Reserved



Independent Medical Associates

www.I-MA.com · info@I-MA.com
p: 888.548.4462 · f: 888-548-1462