

**TARGETED DELIVERY OF
MEDICATIONS AND FLUIDS**

EBI
A BIOMET COMPANY

**The Best Combination of
Durability,**



Steerability,



and Versatility.



EBI®
TargetCath™
System

Fluoro-Guided Steerable Catheter System



Fluoro-Guided epidural pain management for chronic low back pain

Targeted Delivery of Approved Medications in the Epidural Space

EBI[®] TargetCath[™] provides a more durable catheter shaft for better control. The enhanced steering mechanism allows a 135° turning range, left to right. This versatile, dual-port system can be used for the delivery of drugs and contrast, as well as for assisting in diagnosis and treatment.

- Target delivery of approved drugs for chronic low back pain under fluoroscopic guidance
- Enhanced, more controllable and further-reaching steering mechanism compared to existing systems
- More durable catheter shaft allows for easy entry and better control while maneuvering and re-directing in the epidural space of the lumbar spine
- A highly effective tool for mechanical blunt dissection in the epidural space
- Dual lumen catheter offers the ability to efficiently connect and deliver both approved drugs and the administration of fluids indicated for the epidural space
- Single-use sterile packaged steerable catheter and access components eliminate cross contamination concerns
- Ergonomically-designed, steerable catheter handle allows easy grip and precise tip placement

Ordering Information

EBI[®] TargetCath[™] Fluoro-Guided Steerable Catheter

Ref.: TC010

Box of five (5) each individually sterile packaged steerable catheters:

- 2.4mm dia. x 28cm radiopaque to the steerable tip
- Dual lumen with Luer Lock port fillings
- Without detents
- Gamma radiation sterilized

EBI[®] TargetCath[™] Fluoro-Guided Steerable Catheter, Disposable Over-the-Needle Accessory Kit

Ref.: TC002

Box of five (5) each individually sterile packaged accessory kits:

- 1- 8.0F x 6cm Introducer Sheath with Dilator and 1mm Tuohy Tip Stylet

EBI[®] TargetCath[™] Fluoro-Guided Steerable Catheter, Disposable Accessory Kit

Ref.: TC004

Box of five (5) each individually sterile packaged accessory kits:

- 1- 17G x 3 1/2 in. Tuohy Spinal Needle
- 1- .032 in. x 92cm Guidewire, Straight Coated
- 1- 8.0F x 6cm Introducer Sheath with Dilator
- 1- 12F x 10cm Introducer Dilator
- Introducer sheaths and dilator radiopaque
- Ethyl Oxide gas sterilized
- Luer Lock introducer port filling

Over-the-Needle Introducer



For more information, please contact us:
www.I-MA.com
info@I-MA.com
Phone: 888.548.4462
Fax: 888.548.1462