





# Silicone Enteral Feeding Tube




*The hub is clearly marked with "ORAL" and French size to readily identify the tube as an enteral feeding tube.*



*The open end design, with precision drilled side holes, enhances fluid flow and decreases potential feeding stagnation and mucosal irritation.*

## Silicone Enteral Feeding Tube

- Made of soft, biocompatible medical grade silicone designed to enhance extended dwell protocols
- Labeling supports ISMP recommendations that all distal ports and tubing of access lines clearly indicate what is the use of the port/line
- Bold, indelible markings every centimeter
- Color coded for patient safety
- Enteral connector designed to prevent misconnections
- Integrated closure plug for easy cleaning
- Radiopaque stripe for confirmation of tip placement by x-ray
- Latex and DEHP free



Ordering Information for:

# Silicone Feeding Tubes

with Enteral Connector



---

---

## Silicone 40cm Enteral Feeding Tubes

S1FTS5.0-EO	Silicone Feeding Tubes	5.0 Fr.	10/cs
S1FTS6.5-EO	Silicone Feeding Tubes	6.5 Fr.	10/cs
S1FTS8.0-EO	Silicone Feeding Tubes	8.0 Fr.	10/cs

---

---

## Silicone 60cm Enteral Feeding Tubes

S1FTM5.0-EO	Silicone Feeding Tubes	5.0 Fr.	10/cs
S1FTM6.5-EO	Silicone Feeding Tubes	6.5 Fr.	10/cs
S1FTM8.0-EO	Silicone Feeding Tubes	8.0 Fr.	10/cs

---

---

## Silicone 90cm Enteral Feeding Tubes

S1FTL5.0-EO	Silicone Feeding Tubes	5.0 Fr.	10/cs
S1FTL6.5-EO	Silicone Feeding Tubes	6.5 Fr.	10/cs
S1FTL8.0-EO	Silicone Feeding Tubes	8.0 Fr.	10/cs

*"Helping tiny babies make their own footprints in the world..."*

Distributed By:



www.I-MA.com  
info@I-MA.com  
p: 888.548.4462  
f: 888-548-1462



1203 Camden  
San Antonio, Texas 78215  
210.226.2600  
www.footprintmed.com