

ABS SYSTEM

ASPIRATION SYRINGE WITH WINGS

- The syringe called **ABS SYSTEM "with stopper"** allows to produce vacuum making easier the biopsy according to the Menghini procedure.
- The **ABS SYRINGE "without stopper"** makes it easier the aspiration manoeuvre during a cytological sampling collection.

Syringe capacity: 12ml

The wings placed on the plunger make unnecessary the use of "gun" devices, allowing a better and safer manoeuvrability during the procedure.

The reduced distance between the hand that controls the syringe and the needle point increases the user's control.

Male luer lock hub that allows to connect needles, extension tubes and other devices.



CE 0373

CODE	Description	Box
ABS20-B	Winged syringe with stopper	25pcs.
ABS20	Winged syringe without stopper	25pcs.