

ALM SET

BREAST LOCALIZATION NEEDLE

With hook "V" point

The **ALM** is a needle for pre-surgery localization of breast non-palpable lesions.

Features:

- Marked needle with a hook "V" point metal wire
- The funnel shaped clearThe original design of the needle point allows to obtain a large sample perfectly cylindrical.
- The stylet with its pyramidal point allows an easy introduction.The depth marks and stopper permit a safe needle positioning.
- The echogenic marker prevents the spreading of cancerous cells.
- The external cannula can be removed.

A.L.M. Set "V" type

CODE	Diameter	Length	Box
ALM2005	20G	50mm	10
ALM20075	20G	75mm	10
ALM2010	20G	100mm	10
ALM2015	20G	150mm	10
ALM2105	21G	50mm	10
ALM21075	21G	75mm	10
ALM2110	21G	100mm	10

Stiffening cannula

CODE		Box
STIF 50	Stiffening cannula for 50mm guide-wire	10
STIF 75	Stiffening cannula for 75mm guide-wire	10
STIF 100	Stiffening cannula for 100mm guide-wire	10

CE0373



Guide-wire with hook "V" point

ALMS

BREAST LOCALIZATION NEEDLE

The **ALMS** is a needle for pre-surgery localization of breast non-palpable lesions.

Features:

- Needle with a hook point ("V") metal wire at the distal tip
- The reduced dimensions of the *SLIMLINE*™ needle's hub allow to get through the holes of the compression plate
- The wire has two depth marks: while the former indicates that the wire tip is placed nearby the needle tip, the latter shows the complete protrusion of the wire tip from the needle
- The stiffening cannula allows to protect the wire during the surgical operation

CODE	Diameter	Length	Box
ALMS2005	20G	50mm	10
ALMS20075	20G	75mm	10
ALMS2010	20G	100mm	10
ALMS2105	21G	50mm	10
ALMS21075	21G	75mm	10
ALMS2110	21G	100mm	10

Stiffening cannula

CODE		Box
STIF 50	Stiffening cannula for 50mm guide-wire	10
STIF 75	Stiffening cannula for 75mm guide-wire	10
STIF 100	Stiffening cannula for 100mm guide-wire	10

CE0373



Guide-wire with hook "V" point

ALM SET

BREAST LOCALIZATION NEEDLE

With "Z" point

The **ALM** is a needle for pre-surgery localization of breast non-palpable lesions.

Features:

- Marked needle with a "Z" point metal wire.
- The funnel shaped clear luer-hub allows an easier wire penetration.
- The wire has two depth marks: while the former indicates that the wire tip is placed nearby the needle tip, the latter shows the complete protrusion of the wire tip from the needle
- The needle allows both coloured liquids injection and fluids aspiration.
- The stiffening cannula allows to protect the wire during the surgical operation.

A.L.M. Set "Z" type

CODE	Diameter	Length	Box
ALMZ2005	20G	50mm	10
ALMZ20075	20G	75mm	10
ALMZ2010	20G	100mm	10
ALMZ2015	20G	150mm	10

Stiffening cannula

CODE		Box
STIF 50	Stiffening cannula for 50mm guide-wire	10
STIF 75	Stiffening cannula for 75mm guide-wire	10
STIF 100	Stiffening cannula for 100mm guide-wire	10

CE0373



Guide-wire with "Z" point

I.M.A. Hospital Services S.p.A. - Via S. Antonio 150 - 40131 - BO - Italy