

BMN "I" TYPE

BONE MARROW ASPIRATION SET

Bone marrow aspiration set for sternal-iliac crest

Features:

- The needle and stylet points are designed to easily penetrate into the bone marrow cavity. The adjustable and removable needle stop, allowing the needle length variation through a micrometric scale, safely controls depth of penetration.
- The ergonomically designed handle enables a safer manoeuvrability, since the strength needed to penetrate the bone marrow cavity is homogeneously distributed on the whole palm of your hand and not locally.

CODE	Diameter	Length	Box
I 1407M	14G	9-20-53-70mm	10
I 1408M	14G	19-30-63-80mm	10
I 1507M	15G	9-20-53-70mm	10
I 1508M	15G	19-30-63-80mm	10
I 1510M	15G	39-50-83-100mm	10
I 1607M	16G	9-20-53-70mm	10
I 1610M	16G	39-50-83-100mm	10
I 1807M	18G	9-20-53-70mm	10
I 1810M	16G	39-50-83-100mm	10



CE0373

BONE MARROW BIOPSY

IMA
Independent Medical Associates

www.I-MA.com
info@I-MA.com
p: 888.548.4462
f: 888-548-1462