

# BIOMOL

MENGHINI TYPE SEMI-AUTOMATIC BIOPSY DEVICE

The **BIOMOL** makes the biopsy procedure automatic, allowing the physician the sampling control.

The spring not only simplifies the biopsy procedure using only one hand, but allows a much faster performance.

The original design of the needle point allows to obtain a large sample perfectly cylindrical.

The stylet with its pyramidal point allows an easy introduction. The depth marks and stopper permit a safe needle positioning.

The echogenic marker prevents the spreading of cancerous cells.

The external cannula can be removed.



CE 0373

CODE	Diameter	Length	Box
BM 1607	16G	70mm	10
BM 1609	16G	90mm	10
BM 1612	16G	120mm	10
BM 1615	16G	150mm	10
BM 1630	16G	300mm	10
BM 1709	17G	90mm	10
BM 1715	17G	150mm	10
BM 1720	17G	200mm	10
BM 1807	18G	70mm	10
BM 1810	18G	100mm	10
BM 1815	18G	150mm	10
BM 1820	18G	200mm	10
BM 1910	19G	100mm	10
BM 1915	19G	150mm	10
BM 1920	19G	200mm	10
BM 2010	20G	100mm	5
BM 2015	20G	150mm	5
BM 2020	20G	200mm	5
BM 2110	21G	100mm	5
BM 2115	21G	150mm	5
BM 2120	21G	200mm	5
BM 2123	21G	230mm	5
BM 2210	22G	100mm	5
BM 2215	22G	150mm	5
BM 2220	22G	200mm	5

**IMA**  
Independent Medical Associates

www.I-MA.com  
info@I-MA.com  
p: 888.548.4462  
f: 888-548-1462

HISTOLOGICAL BIOPSY