

ECOJEKT

ASPIRATION NEEDLE WITH MODIFIED CHIBA POINT

The **ECOJEKT** is a needle with a modified Chiba type point (double sharpness).

This needle includes a removable stylet perfectly coupled to the needle.

Both the needle sliding stopper and the depth marks allow to safely control the penetration depth. The inner echogenic marker prevents the spreading of cancerous cells.

Clear luer hub with colour code.

The **ECOJEKT** is a universal needle. Its applications are: biopsies, aspirations, injections of contrast media, amniotic fluid sampling, alcoholizations, exploratory punctures, withdrawal of chorial villi, lumbar punctures.

H.S. Hospital Services S.p.A. is certified ISO 9001 and EN 45001

CODE	Diameter	Length	Box	CODE	Diameter	Length	Box
E1715	17G	150mm	10	E2035	20G	350mm	10
E1720	17G	200mm	10	E2107	21G	70mm	10
E1815	18G	150mm	10	E2110	21G	100mm	10
E1818	18G	180mm	10	E2115	21G	150mm	10
E1820	18G	200mm	10	E2118	21G	180mm	10
E1825	18G	250mm	10	E2120	21G	200mm	10
E1828	18G	280mm	10	E2130	21G	300mm	10
E1907	19G	70mm	10	E2205	22G	50mm	10
E1910	19G	100mm	10	E22075	22G	75mm	10
E1915	19G	150mm	10	E2210	22G	100mm	10
E1918	19G	180mm	10	E2215	22G	150mm	10
E1920	19G	200mm	10	E2218	22G	180mm	10
E1930	19G	300mm	10	E2220	22G	200mm	10
E2005	20G	50mm	10	E2315	23G	150mm	10
E2007	20G	70mm	10	E2318	23G	180mm	10
E2010	20G	100mm	10	E2320	23G	200mm	10
E2015	20G	150mm	10	E2507*	25G	70mm	10
E2018	20G	180mm	10	E2510*	25G	100mm	10
E2020	20G	200mm	10	E2515*	25G	150mm	10
E2030	20G	300mm	10				

*Codes with external echogenic marker

CE 0373

IMA

Independent Medical Associates

www.I-MA.com

info@I-MA.com

p: 888.548.4462

f: 888-548-1462

CYTOLOGICAL BIOPSY