

# HS BIOPSY PUNCH

S K I N B I O P S Y D E V I C E

Using the **BIOPSY PUNCH**, surgical removal of skin small cancerous cells has become easier, faster, safer and effective.

#### BIOPSY PUNCH - Skin Biopsy

- Tissue sampling
- Check of the exeresis integrity
- Determination of the RSTL (relaxed Skin tension Lines)

CODE	Diameter (mm)	Box
HBP 20	2	20
HBP 25	2.5	20
HBP 30	3	20
HBP 35	3.5	20
HBP 40	4	20
HBP 50	5	20
HBP 60	6	20
HBP 80	8	20

- Colour code
- Size of the diameter impressed on the hub
- Point protective cap



CE 0373



HISTOLOGICAL BIOPSIY

H.S. Hospital Services S.p.A. is certified ISO 9001 and EN 46001