

HS TRAPSYSTEM® SET

BONE MARROW BIOPSY SET WITH TRAPSYSTEM DEVICE

The **TRAPSYSTEM®** grants a quality marrow sampling by trapping the specimen into the needle cannula without damaging it and without performing the painful needle deflection technique.

- Greater patient comfort
- Safer marrow sampling
- Simpler manoeuvrability
- Faster procedure

The needle, with its concave cutting edge point, allows to easily penetrate into the bone marrow cavity. Also the marked stylet, with its pyramidal point, allows to easily penetrate into the bone marrow cavity. The ergonomically designed handle enables a safer manoeuvrability, since the strength needed to penetrate the bone marrow cavity is homogeneously distributed on the whole palm of your hand and not locally.

CODE	Diameter	Length	Box
TRAPJ 0810	8G	100mm	10
TRAPJ 0815	8G	150mm	10
TRAPJ 1110	11G	100mm	10
TRAPJ 1115	11G	150mm	10
TRAPJ 1306	13G	65mm	10
TRAPJ1310	13G	100mm	10



HS Medical Services S.p.A. - certified ISO 9001 and ISO 13485

IMA
Independent Medical Associates

11733 66th St N Ste 113
Largo, FL 33773
Ph: 888-548-4462
Fax: 888-548-1462
info@I-MA.com
www.I-MA.com

BONE MARROW BIOPSY