

MB SET

DEVICE FOR BREAST AND THYROID BIOPSY

"Lancet" pointed needle that easily penetrates the tissue without lacerating it.

The clear luer hub allows the visual control both of the sampling collection and of its macroscopic features.

The stopper makes it possible to safely control the penetration depth.

The ABS type syringe makes unnecessary the use of "gun" devices offering anyway, with its wings, a better and safer manoeuvrability.

The reduced distance between the hand that manoeuvres the syringe and the needle's point increases both sensitivity and precision in centring the nodule.

CODE	Diameter	Diameter	Length	Box
MB 2005	20G	0.88mm	50mm	20
MB 2007	20G	0.88mm	70mm	20
MB 2205	22G	0.71mm	50mm	20
MB 2207	22G	0.71mm	70mm	20
MB 2505	25G	0.51mm	50mm	20



CE 0373

H.S. Hospital/Surviva S.p.A. a certified ISO 9001 and EN 14961

IMA
Independent Medical Associates

www.I-MA.com
info@I-MA.com
p: 888.548.4462
f: 888-548-1462

CYTOLOGICAL BIOPSY