

SECURCUT

MENGHINI TYPE MANUAL BIOPSY DEVICE

The **SECURCUT** is a biopsy needle for soft tissue. The sampling technique is the classic modified Menghini type.

The original design of the needle point allows to obtain a large sample perfectly cylindrical.

The stylet, with its pyramidal tip coupled with the syringe plunger, allows an easy penetration.

The depth marks and stopper permit a safe positioning of the needle.

The internal echogenic marker prevents the spreading of cancerous cells.

The external cannula can be removed.

H.S. Hospital Service S.p.A. is certified ISO 9001 and EN 45001

CODE	Diameter	Length	Box	CODE	Diameter	Length	Box
S2207	22G	70mm	10	S1807	18G	70mm	10
S2210	22G	100mm	10	S1809	18G	90mm	10
S2212	22G	120mm	10	S1810	18G	100mm	10
S2215	22G	150mm	10	S1812	18G	120mm	10
S2220	22G	200mm	10	S1815	18G	150mm	10
S2104	21G	40mm	10	S1820	18G	200mm	10
S2107	21G	70mm	10	S1707	17G	70mm	10
S2109	21G	90mm	10	S1709	17G	90mm	10
S2110	21G	100mm	10	S1710	17G	100mm	10
S2112	21G	120mm	10	S1712	17G	120mm	10
S2115	21G	150mm	10	S1715	17G	150mm	10
S2120	21G	200mm	10	S1720	17G	200mm	10
S2009	20G	90mm	10	S1607	16G	70mm	10
S2012	20G	120mm	10	S1609	16G	90mm	10
S2015	20G	150mm	10	S1610	16G	100mm	10
S2020	20G	200mm	10	S1612	16G	120mm	10
S1905	19G	50mm	10	S1615	16G	150mm	10
S1907	19G	70mm	10	S1620	16G	200mm	10
S1909	19G	90mm	10	S1509	15G	90mm	10
S1910	19G	100mm	10	S1510	15G	100mm	10
S1915	19G	150mm	10	S1512	15G	120mm	10
S1920	19G	200mm	10	S1515	15G	150mm	10

CE 0373

IMA
Independent Medical Associates

www.I-MA.com
info@I-MA.com
p: 888.548.4462
f: 888-548-1462



H.S. HOSPITAL SERVICE S.p.A.
Via delle Valli snc - 04011 Aprilia (LT)
Tel. 089201901 - Fax 0892727871 - Orders fax 089276519
www.hs-hospital-service.com - e-mail: hs@hs-hospital-service.com

HISTOLOGICAL BIOPSY