

TOP CUT M

TRU-CUT NEEDLE FOR AUTOMATIC DEVICES

The **TOP CUT** is a tru-cut biopsy needle for **MAGNUM** type automatic systems.

The **TOP CUT** needles features are:

- Marked needle with Menghini type tru-cut point
- Whistle pointed stylet
- Inner echogenic marker
- 16cm notch
- 22mm advancement
- Colour code

TOP CUT - 22mm advancement

CODE	Diameter	Length	Box
TOP 1410 M	14G	100mm	10
TOP 1413 M	14G	130mm	10
TOP 1416 M	14G	160mm	10
TOP 1420 M	14G	200mm	10
TOP 1610 M	16G	100mm	10
TOP 1616 M	16G	160mm	10
TOP 1620 M	16G	200mm	10
TOP 1810 M	18G	100mm	10
TOP 1816 M	18G	160mm	10
TOP 1820 M	18G	200mm	10
TOP 1825 M	18G	250mm	10
TOP 1830 M	18G	300mm	10
TOP 2010 M	20G	100mm	10
TOP 2016 M	20G	160mm	10
TOP 2020 M	20G	200mm	10

*MAGNUM is a Bard INC.USA registered trade-mark



CE 0373

HISTOLOGICAL BIOPSIY

IMA
Independent Medical Associates

www.I-MA.com
info@I-MA.com
p: 888.548.4462
f: 888-548-1462

HS

H.S. HOSPITAL SERVICE S.p.A.
Via delle Valli snc - 04011 Aprilia (LT)
Tel. 069201961 - Fax 0692727871 - Orders fax 069275519
www.hshospitalservice.com e-mail: hs@hshospitalservice.com