





## CATALOG



## **ACCU-DISC™**

**Fluid Dispensing Syringe**

- The ACCU-DISC™ Pressure Monitoring System offers the advantage of accurate delivery of fluids to the body with the ability to monitor the pressure of those fluids in discography interpretation.



# ANCILLARY ITEMS

Stock #	Description
<b>Quincke Point Spinal Needles</b>	
PISN1835	18G x 3.5" Quincke
PISN2035	20G x 3.5" Quincke
PISN2235	22G x 3.5" Quincke
PISN2335	23G x 3.5" Quincke
PISN2525	25G x 2.5" Quincke
PISN2535	25G x 3.5" Quincke
PISN2735	27G x 3.5" Quincke

<b>Quincke Point Long Spinal Needles</b>	
PISN2050	20G x 5" Quincke
PISN2060	20G x 6" Quincke
PISN2250	22G x 5" Quincke
PISN2260	22G x 6" Quincke
PISN2270	22G x 7" Quincke
PISN2280	22G x 8" Quincke
PISN2350	23G x 5" Quincke
PISN2360	23G x 6" Quincke
PISN2550	25G x 5" Quincke
PISN2560	25G x 6" Quincke
PISN2570	25G x 7" Quincke
PISN2580	25G x 8" Quincke

Stock #	Description
<b>Pencil Point Spinal Needles</b>	
PIPP2235	22G x 3.5"
PIPP2535	25G x 3.5"

<b>Chiba Needles</b>	
PICH2235	22G x 3.5"
PICH2060	20G x 6"
PICH2250	22G x 5"
PICH2260	22G x 6"
PICH2270	22G X 7"
PICH2280	22G X 8"

<b>Curved Chiba Needles</b>		
2260 CCT	Curved Tip Chiba Needle	22G x 6"
8600 CT	Curved Tip Chiba Needle w/ Introducer Needle	22G x 6" w/ 18G x 6.5" Introducer
2235BC	Blunt Tipped Curved w/ Introducer Needle	22G x 3.5" w/ 16G x 2.0" Introducer
2260BC	Blunt Tipped Curved w/ Introducer Needle	22G x 6.0" w/ 16G x 3.25" Introducer

<b>Epidural Steroid Injection Needles</b>	
2025T	20G X 2.5"
PITRW1735	17G x 3.5"
PITRW1835	18G x 3.5"
PITRW2035	20G x 3.5"
PITRW2235	22G x 3.5"



# ANCILLARY ITEMS

Stock #	Description
<b>Long Epidural Steroid Injection Needles</b>	

PITRW1750	17G x 5"
1760T	17G x 6"
PITRW1850	18G x 5"
PITRW2050	20G x 5"

<b>Loss of Resistance Syringes</b>	
PISLL7	Plastic LOR Syringe 7CC LL
PISLS7	Plastic LOR Syringe 7CC LS
PISGLL5	Glass LOR Syringe 5CC LL
PISGLS5	Glass LOR Syringe 5CC LS

PISLL7	Plastic LOR Syringe 7CC LL
PISLS7	Plastic LOR Syringe 7CC LS
PISGLL5	Glass LOR Syringe 5CC LL
PISGLS5	Glass LOR Syringe 5CC LS

<b>Catheters</b>	
CA20CES	20G Close End Catheter w/Stylet
CA20OE	20G Open End Catheter
CA20CE	20G Close End Catheter

CA20CES	20G Close End Catheter w/Stylet
CA20OE	20G Open End Catheter
CA20CE	20G Close End Catheter

<b>Extension Sets</b>	
CA12T	Bi-fercated extension set 12"
CUS1065	6" Mini-Bore Extension Set
CUS1027	12" Mini - Bore Extension Set

CA12T	Bi-fercated extension set 12"
CUS1065	6" Mini-Bore Extension Set
CUS1027	12" Mini - Bore Extension Set

Stock #	Description	Size
<b>Radiation Reduction Gloves</b>		

RR1 - 6.5	6.5
RR1 - 7.0	7.0
RR1 - 7.5	7.5
RR1 - 8.0	8.0
RR1 - 8.5	8.5
RR1 - 9.0	9.0
RR2 - 7.0	7.0
RR2 - 7.5	7.5
RR2 - 8.0	8.0
RR2 - 8.5	8.5
RR2 - 9.0	9.0

RR1 60KVP 45% (.008), RR2 60KVP 55% (.008), THINNER FINGER TIPS  
ANATONICALLY CORRECT CURVED FINGERS, OTHER SIZES AVAILABLE

<b>Discogram Needles</b>	
DDN5000	18G x 6 3/16" Cannula
	20G x 6 1/8" sharp trocar needle
	20G x 6 1/8" blunt obturator needle
	22 Gauge injection needle
8620	18G x 6" Spinal needle w/ 22G x 8" Spinal needle
8630	18G x 6-3/16" Spinal needle w/ 20G x 6-1/8" Spinal needle
8600	18G x 3.5" Spinal needle w/ 22G x 6" Spinal needle
8610	20G x 3.5" Spinal needle w/ 25G x 6" Spinal needle

DDN5000	18G x 6 3/16" Cannula
	20G x 6 1/8" sharp trocar needle
	20G x 6 1/8" blunt obturator needle
	22 Gauge injection needle
8620	18G x 6" Spinal needle w/ 22G x 8" Spinal needle
8630	18G x 6-3/16" Spinal needle w/ 20G x 6-1/8" Spinal needle
8600	18G x 3.5" Spinal needle w/ 22G x 6" Spinal needle
8610	20G x 3.5" Spinal needle w/ 25G x 6" Spinal needle



# ANCILLARY ITEMS

Stock #	Description
<b>Radio Frequency Needles</b>	
PIRF20-100-10CS	RF CANNULAE 20G X 100MM 10MM ACTIVE TIP CURVED SHARP
PIRF20-150-10CS	RF CANNULAE 20G X 150MM 10MM ACTIVE TIP CURVED SHARP
PIRF20-100-10CB	RF CANNULAE 20G X 100MM 10MM ACTIVE TIP CURVED BLUNT
PIRF20-150-10CB	RF CANNULAE 20G X 150MM 10MM ACTIVE TIP CURVED BLUNT
PIRF22-54-4	RF CANNULAE 22G X 54MM 4MM ACTIVE TIP STRAIGHT
PIRF22-100-5	RF CANNULAE 22G X 100MM 5MM ACTIVE TIP STRAIGHT
PIRF20-145-5	RF CANNULAE 20G X 145MM 5MM ACTIVE TIP STRAIGHT
PIRF21-100-5	RF CANNULAE 21G X 100MM 5MM ACTIVE TIP STRAIGHT
PIRF22-100-10	RF CANNULAE 22G X 100MM 10MM ACTIVE TIP STRAIGHT
PIRF20-145-10	RF CANNULAE 20G X 145MM 10MM ACTIVE TIP STRAIGHT
PIRF20-145-15	RF CANNULAE 20G X 145MM 15MM ACTIVE TIP STRAIGHT
PIRF23-100-5	RF CANNULAE 23G X 100MM 5MM ACTIVE TIP STRAIGHT
NP20145T	RF CANNULAE 20G X 145MM 10MM ACTIVE TIP STRAIGHT
NP2210010	RF CANNULAE 22G X 100MM 10MM ACTIVE TIP STRAIGHT
NP22100	RF CANNULAE 22G X 100MM 5MM ACTIVE TIP STRAIGHT
PIRF22-100-10CS	RF CANNULAE 22G X 100MM 10MM ACTIVE TIP CURVED SHARP
PIRF18-100-10CS	RF CANNULAE 18G X 100MM 10MM ACTIVE TIP CURVED SHARP
PIRF18-100-5CS	RF CANNULAE 18G X 100MM 5MM ACTIVE TIP CURVED SHARP
PIRF18-150-10CS	RF CANNULAE 18G X 150MM 10MM ACTIVE TIP CURVED SHARP

# STOCK TRAY

## Universal Block Tray

### Catalog Number CUS803

#### Qty. Description

- |   |                          |
|---|--------------------------|
| 2 | Syringe 10CC (Luer Lock) |
| 1 | Needle 21G X 1-1/2"      |
| 1 | Needle 25G X 1-1/2"      |
| 3 | Gauze Sponge 4 X 4-8 Ply |
| 3 | Sponge Applicators       |
| 1 | Povidone Iodine 1 oz.    |
| 1 | Fenestrated Drape        |
| 1 | Hospital Wrap            |
| 1 | Tray, Large Deep         |

### Catalog Number CUS1428

#### Qty. Description

- |   |   |
|---|---|
| 1 | Needle 18G X 1-1/2"                               |
| 1 | Syringe 5CC (Luer Lock)                           |
| 1 | Syringe 10CC (Luer Lock)                          |
| 1 | Syringe 3CC (Luer Lock)                           |
| 1 | Needle 25G X 1-1/2"                               |
| 1 | Sodium Chloride 0.9% 10ML Ampule                  |
| 1 | Lidocaine 1% 5ML Ampule                           |
| 1 | Povidone Iodine 1 oz.                             |
| 3 | Sponge Applicators                                |
| 1 | Absorbent Towel                                   |
| 1 | Clear Fenestrated Drape, Tape Strips              |
| 4 | Gauze Sponge 4 X 4-12 Ply                         |
| 1 | Hospital Wrap                                     |
| 1 | Extension Set 12" Smallbore, Clamp, Rotating Luer |

### Catalog Number CUS610

#### Qty. Description

- |   |                          |
|---|--------------------------|
| 1 | Syringe 10CC (Luer Lock) |
| 2 | Syringe 5CC (Luer Lock)  |
| 1 | Syringe 3CC (Luer Lock)  |
| 1 | Needle 18G X 1-1/2"      |
| 1 | Needle 27G X 1-1/4"      |
| 3 | Sponge Applicators       |
| 3 | Gauze Sponge 4 X 4-8 Ply |
| 1 | Povidone Iodine 1 oz.    |
| 1 | Lidocaine 1% 5ML Ampule  |
| 1 | Fenestrated Drape        |
| 1 | Needle Stick Pad         |
| 1 | Hospital Wrap            |

# STOCK TRAY

## Single Shot Trays, Glass LOR

### Catalog Number CUS001

#### Qty. Description

1	Tuohy Needle 18G X 3-1/2"
1	Needle 22G X 1-1/2"
1	Needle 18G X 1-1/2"
1	Needle 27G X 1-1/4"
1	Filter Needle 19G X 1 1/2"
1	Syringe 10CC (Luer Lock)
1	Syringe 3CC (Luer Lock)
1	Glass Syringe 5CC (Luer Lock)
1	Sodium Chloride .9% 10CC Ampule
1	Lidocaine 1% 5ML Ampule
4	4 X 4-8 Ply Gauze
1	Absorbent Towel
1	Fenestrated Drape W/Tape
1	Povidone Iodine Triple Swabs
1	Large Alcohol Wipe
1	Hospital Wrap

### Catalog Number CUS002

#### Qty. Description

1	Tuohy Needle 20G X 3-1/2"
1	Needle 22G X 1-1/2"
1	Needle 18G X 1-1/2"
1	Needle 27G X 1-1/4"
1	Filter Needle 19G X 1 1/2"
1	Syringe 10CC (Luer Lock)
1	Syringe 3CC (Luer Lock)
1	Glass Syringe 5CC (Luer Lock)
1	Sodium Chloride .9% 10CC Ampule
1	Lidocaine 1% 5ML Ampule
4	4 X 4-8 Ply Gauze
1	Absorbent Towel
1	Fenestrated Drape W/Tape
1	Povidone Iodine Triple Swabs
1	Large Alcohol Wipe
1	Hospital Wrap

### Catalog Number CUS1517

#### Qty. Description

1	Needle 18G X 1-1/2"
1	Syringe 5CC (Luer Lock)
1	Syringe 10CC (Luer Lock)
1	Glass Syringe 5CC, Metal (Luer Lock)
1	Syringe 3CC (Luer Lock)
1	Needle 25G X 1-1/2"
1	Tuohy Epidural Needle 20G X 3-1/2", Fixed Wing, Metal Stylet
1	Sodium Chloride 0.9% 10ML Ampule
1	Lidocaine 1% 5ML Ampule
1	Povidone Iodine 1oz.
3	Sponge Applicators
1	Absorbent Towel
1	Clear Fenestrated Drape, Tape Strips
4	Gauze Sponge 4 X 4-12 Ply
1	Hospital Wrap

### Catalog Number CUS1427

#### Qty. Description

1	Needle 18G X 1-1/2"
1	Syringe 5CC (Luer Lock)
1	Syringe 10CC (Luer Lock)
1	Glass Syringe 5CC, Metal (Luer Lock)
1	Syringe 3CC (Luer Lock)
1	Needle 25G X 1-1/2"
1	Tuohy Epidural Needle, 18G X 3-1/2" R/W Metal Stylet
1	Sodium Chloride 0.9% 10ML Ampule
1	Lidocaine 1% 5ML Ampule
1	Povidone Iodine 1oz.
3	Sponge Applicators
1	Absorbent Towel
1	Clear Fenestrated Drape, Tape Strips
4	Gauze Sponge 4 X 4-12 Ply
1	Hospital Wrap
1	Header Bag 10 X 15



# STOCK TRAY

## Single Shot Trays, Plastic LOR

### Catalog Number CUS096

#### Qty. Description

1	Needle 18G x 1.50"
1	Needle 27G x 1.25"
1	Syringe 5CC (Luer lock)
1	Syringe 10CC (Luer lock)
1	Plastic L.O.R. Syringe 7CC (Luer Slip)
1	Tuohy Epidural Needle 20G x 3.50", Wings
1	Needle 21G x 1.50"
1	Sodium chloride 0.9% 10ML ampule
2	Lidocaine 1% 5ML Ampule
1	Povidone Iodine 1oz.
3	Sponge Applicators
3	Gauze Sponge 4 x 4-8 Ply
1	Absorbent Towel
1	Fenestrated Drape
1	Two Deck Pain Tray
1	Hospital Wrap

### Catalog Number CUS1243

#### Qty. Description

1	Tuohy Epidural Needle 22G X 3-1/2"
1	Syringe 5CC (Luer Lock)
1	Syringe 10CC (Luer Lock)
2	Syringe 3CC (Luer Lock)
1	Needle 18G X 1-1/2"
1	Needle 25G X 1-1/2"
1	Plastic L.O.R. Syringe 7CC (Luer Slip)
1	Sodium Chloride 0.9% 10ML Ampule
1	Lidocaine 1% 5ML Ampule
4	Gauze Sponge 4 X 4-8 Ply
1	Foam Syringe Pouch
1	Opaque Drape, 4inch Fenestration 24 X 24, Tape Strips
1	Povidone Iodine 1oz.
3	Sponge Applicators
2	Graduated Medicine Cup 2oz.
1	Two Compartment Tray
1	Hospital Wrap

### Catalog Number CUS628

#### Qty. Description

2	Needle 18G x 1.50"
1	Needle 27G x 1.25"
1	Syringe 5CC (Luer Lock)
1	Syringe 10cc (Luer Lock)
1	Syringe 3cc (Luer Lock)
1	Tuohy Epidural Needle 18G x 3.50"
1	Plastic L.O.R. Syringe 7CC (Luer Lock)
1	Sodium Chloride 0.9% 10ML Ampule
2	Lidocaine 1% 5ML Ampule
4	Gauze Sponge 4 x 4-8 ply
1	Foam Syringe Pouch
1	Fenestrated Drape
1	Graduated Medicine Cup 2oz.
1	Large Deep Prep Tray 9 x 5 x 2"
1	Hospital Wrap

### Catalog Number CUS025

#### Qty. Description

1	Needle 18G X 1-1/2"
1	Syringe 5CC (Luer Lock)
1	Syringe 10CC (Luer Lock)
1	Syringe 3CC (Luer Lock)
1	Needle 25G X 1-1/2"
1	Plastic L.O.R. Syringe 7CC (Luer Lock)
1	Tuohy Epidural Needle, 18G X 3-1/2" R/W Metal Stylet
1	Sodium Chloride 0.9% 10ML Ampule
1	Lidocaine 1% 5ML Ampule
1	Povidone Iodine 1oz.
3	Sponge Applicators
4	Gauze Sponge 4 X 4-8 Ply
1	Fenestrated Drape
1	Hospital Wrap
1	Header Bag 10 X 15

# STOCK TRAY

## Nerve Block Tray

### Catalog Number CUS626

Qty.	Description
1	Needle 18G X 1.50"
1	Syringe 5CC (Luer Lock)
1	Syringe 20CC (Luer Lock)
1	Needle 25G X 1.50"
1	Syringe 1CC (Luer Slip)
1	Quincke Spinal Needle 22G X 3.50"
1	Sodium Chloride 0.9% 10ML Ampule
1	Lidocaine 1% 5ML Ampule
1	Povidone Iodine 1oz.
3	Sponge Applicators
4	Gauze Sponge 4 X 4-8 Ply
1	Adhesive Bandage
1	Clear Fenestrated Drape, Tape Strips
1	Two Deck Pain Tray
1	Hospital Wrap
1	Extension Set 6" Smallbore, Slide Clamp

## Lumbar Puncture Tray

### Catalog Number CUS4301

Qty.	Description
1	Needle 22G X 1-1/2"
1	Needle 25G X 5/8"
1	Syringe 3CC (Luer Lock)
1	Quincke Spinal Needle 20G X 3-1/2"
1	Lidocaine 1% 2ML Ampule
3	Sponge Applicators
1	Adhesive Bandage
3	Gauze Sponge 2 X 2-12 Ply
1	Absorbent Towel
1	Fenestrated Drape
1	Lumbar Puncture Myelogram Tray
1	Clear Specimen Tube #1
1	Clear Specimen Tube #2
1	Clear Specimen Tube #3
1	Clear Specimen Tube #4
1	Hospital Wrap
1	Header Bag 12 X 22
1	Stopcock 3-Way Low Pressure
1	Extension Tubing 6"
1	Manometer 2-Piece 550MM

## Myelogram Tray

### Catalog Number CUS04035

Qty.	Description		Description
1	Needle 18G X 1-1/2"	1	Myelogram Drape W/ 2-3/4" Taped Fen
1	Syringe 5CC (Luer Lock)	1	Lumbar Puncture Myelogram Tray
1	Needle 22G X 1-1/2"	1	Clear Specimen Tube #1
1	Syringe 20CC (Luer Lock)	1	Clear Specimen Tube #2
1	Needle 25G X 5/8"	1	Clear Specimen Tube #3
1	Quincke Spinal Needle 22G X 3-1/2"	1	Hospital Wrap
1	Lidocaine 1% 5ML Ampule	1	Extension Tubing 6"
3	Sponge Applicators	1	Extension Tubing 20"
1	Adhesive Bandage		
3	Gauze Sponge 3 X 3-12 Ply		
1	Absorbent Towel		

# CREATE-A-KIT

Step 1: Complete the following company information form available at [www.integra-ipm.com](http://www.integra-ipm.com)

Company/Customer Name: \_\_\_\_\_

Address: \_\_\_\_\_

City, State, Zip: \_\_\_\_\_

Telephone: \_\_\_\_\_

Fax: \_\_\_\_\_

Email: \_\_\_\_\_

Step 2: Sample

Non Sterile: \_\_\_\_\_

Sterile: \_\_\_\_\_

Step 3: Indicate the Quantity of Desired Products.

## Epidural Needles

- ☐ 17GA x 3-1/2"
- ☐ 18GA x 3-1/2"
- ☐ 18GA x 6"
- ☐ 20GA x 2-1/2"
- ☐ 20GA x 3-1/2"
- ☐ 22GA x 3-1/2"
- ☐ Removable Wing
- ☐ Fixed Wing
- ☐ Huestad
- ☐ Metal Hub
- ☐ Plastic Stylet
- ☐ Metal Stylet
- ☐ Other \_\_\_\_\_

## Loss of Resistance Syringes

- ☐ Plastic Syringe 7CC, (LL)
- ☐ Plastic Syringe 7CC, (LS)
- ☐ Glass Syringe 5CC, (LL)
- ☐ Glass Syringe 5CC, (LS)
- ☐ Other \_\_\_\_\_

## Syringes

- ☐ Syringe 3CC, (LL)
- ☐ Syringe 3CC, (LS)
- ☐ Syringe 5CC, (LL)
- ☐ Syringe 5CC, (LS)
- ☐ Syringe 10CC, (LL)
- ☐ Syringe 10CC, (LS)
- ☐ Syringe 20CC, (LL)
- ☐ Syringe 20CC, (LS)
- ☐ Syringe 1CC, (LS)
- ☐ Other \_\_\_\_\_

## Needles

- ☐ Needle 18GA x 1-1/2" (Pink)
- ☐ Filter Needle 19GA x 1-1/2" (Clear)
- ☐ Needle 20GA x 1-1/2" (Blue)
- ☐ Needle 21GA x 1-1/2" (Green)
- ☐ Needle 22GA x 1-1/2" (Black)
- ☐ Needle 25GA x 5/8" (Blue)
- ☐ Needle 25GA x 1-1/2" (Blue)
- ☐ Needle 27GA x 1/4" (Gray)
- ☐ Filter Straw
- ☐ Other \_\_\_\_\_

## Spinal Needles

- ☐ 18GA x 3-1/2"
- ☐ 20GA x 3-1/2"
- ☐ 22GA x 3-1/2"
- ☐ 22GA x 5"
- ☐ 22GA x 6"
- ☐ 22GA x 7"
- ☐ 22GA x 8"
- ☐ 25GA x 3-1/2"
- ☐ Other \_\_\_\_\_

## Medications

- ☐ Sodium Chloride 0.9% 10ML Ampule
- ☐ Lidocaine 1% 5ML Ampule
- ☐ Lidocaine 1.5% w/1:200. EPI Ampule
- ☐ Lidocaine 1.5% 20ML Ampule
- ☐ Lidocaine 2% 10ML Ampule
- ☐ Bupivacaine .25% 30ML
- ☐ Marcaine 7.5MG w/82.5MG Dextrose 2ML Ampule
- ☐ Other \_\_\_\_\_

## Drapes & Towels

- ☐ Fenestrated Drape
- ☐ Fenestrated Drape w/Tape
- ☐ Utility Drape w/Tape
- ☐ Clear Fenestrated Drape
- ☐ Clear Fenestrated Drape w/Adhesive
- ☐ Absorbent Towel
- ☐ Blue O.R. Towels
- ☐ C-ARM Drape
- ☐ Other \_\_\_\_\_

## Prep Products

- ☐ Sponge Applicators
- ☐ Povidone Iodine 1 oz. Packet
- ☐ Povidone Iodine Triple Swabs
- ☐ Pair Surgeon Gloves
- ☐ Gauze 4" x 4" 8 Ply
- ☐ Gauze 3" x 3" 12 Ply
- ☐ Adhesive Bandage
- ☐ Needle Stick Pad
- ☐ Chloraprep
- ☐ Other \_\_\_\_\_

## Extension Set

- ☐ Extension Set 6" Minibore
- ☐ Extension Set 12" Minibore
- ☐ Extension Set 20" Minibore
- ☐ Other \_\_\_\_\_

## Formed Trays

- ☐ Large Shallow Prep
- ☐ Single Deck Tray
- ☐ Two Deck Pain Tray
- ☐ Other \_\_\_\_\_

## Continuous Epidural Supplies


- ☐ 20GA Epidural Catheter
- ☐ Flat Filter
- ☐ Other \_\_\_\_\_

## Miscellaneous

- ☐ Skin Marker
- ☐ Inflation Syringe

Notes: \_\_\_\_\_







# ACCUMETER

**The Precise Device for Discography**

- Integrated high resolution digital display
- Clear unobstructed view of fluid path
- Optional free plunger movement for tactile feedback
- Smooth continuous plunger advancement
- Ergonomic design for natural hand position





**IMA**  
Independent Medical Associates

11733 66th Street North  
Suite 113 Largo, FL 33773  
Phone: 888-548-4462  
Fax: 888-548-1462  
E-Mail: [Info@I-MA.com](mailto:Info@I-MA.com)  
Web: [www.I-MA.com](http://www.I-MA.com)