





Life-Tech, Inc

# eZono™ 3000 Ultrasound for Anesthesia




- Portable & robust
- B-mode & Color Power Doppler
- Image & Video storage to USB
- Touch screen ease of use



## eZono Cue Cards\*

- Standardized workflow and image interpretation
- Accelerated learning curve on nerve blocks & vascular access
- Enhanced clinical education

# Design and Innovation for Regional Anesthesia



## Image Styles

Exclusive to the eZono 3000, four user selectable image styles




## eZono 3000 Probes

Offering a full range of linear and convex probes to meet the needs of regional anesthesia

**Hard Probe Surface**  
Resistant to needle points, alcohol based disinfectants, adhesive probe covers, and being dropped

**Extra Long Cable**  
The 2 meter cables are ideal for cross bed procedures

**Small Form Factor**  
Simplifies in-plane injections



## eZono 3000 Specifications

### System:

- Weight: 9.7 lb / 4.4 kg
- Dimensions: 12.76" x 11.1" x 3.66" / 32.4 cm x 28.2 cm x 9.3 cm
- 4 GB internal memory
- 10.4" / 26.4 cm diagonal LCD touch screen (24-bit, 60 million colors)
- Digital broadband architecture
- 2 x USB 2.0 ports, 1 Express Card port, built-in microphone

### Software:

- Proprietary operating system
- Less than 20 seconds start time
- Icon-based graphical user interface via touch screen
- Software enhancements available via USB connection throughout the life of the product
- Store up to 30,000 images or 30 minutes of continuous video

Distributed By:

\*Cart optional

**IMA**  
Independent Medical Associates

www.I-MA.com · info@I-MA.com  
p: 888.548.4462 · f: 888-548-1462

