

Infu-Stat™

Disposable Pressure Infuser



- Latex Free
- Loads easily.
- One hand operation for inflation or deflation.
- Large bore stopcock for rapid inflation and deflation.
- Color-coded gauge is easy to read and maintains reliable pressure.
- The safety valve eliminates the possibility of over inflation and releases air at 325mm/hg.
- Sizing is color-coded for easy identification-Graphics are sharp and easy to read.
- Economical, single patient use to prevent cross contamination and eliminates costly cleaning and



Item #	Size	Units
MTM 305	500ml	5/box - 25/cs
MTM 310	1000ml	5/box - 25/cs
MTM 330E	3000ml	2/box - 10/cs

For additional information or to place an order, please feel free to contact:

Matt Scussel

INDEPENDENT MEDICAL ASSOCIATES

11733-66th St. North, Suite 113, Largo, FL 33773

888-548-4462 / 727-548-4462 x734

Fax: 888-548-1462 Cell: 321-662-7287

Email: mscussel@j-ma.com

IMA
Independent Medical Associates



We accept
Visa and Mastercard



Infu-Stat™

Disposable Pressure Infuser

- **100% Latex Free**
- **Cost Competitive**
- **User Friendly**
- **One hand operation**



IMA
Independent Medical Associates

For more information or to place an order, please contact:
Independent Medical Associates
11733 66th Street North
Suite 113
Largo, FL 33773

Phone: 888-548-4462
Fax: 888-548-1462
E-Mail: Info@I-MA.com
Web: www.I-MA.com