



Conventional Laryngoscopes Handles & Locks

Chrome plated and stainless steel waterproof handles are available from SunMed for use with conventional laryngoscope blades. Conventional Howland Lock also available.

For detailed product information, visit www.Sun-Med.com or consult a SunMed product catalog.



Chrome Plated Brass Handles


- Chrome plated brass
- Batteries are not included

Stainless Steel Waterproof Handles

- Stainless steel waterproof handles may be immersed in disinfectant without removing batteries
- Batteries are not included

Cat#	Description	Size	Batteries	Box
5-0237-01	Penlite Handle	143mm	2 "AA"	1
5-0237-02	Stubby Handle	110mm	2 "AA"	1
5-0237-03	Medium Handle	145mm	2 "C"	1
5-0237-05	Large Handle	164mm	2 "D"	1
5-0237-08	Mini Handle	110mm	Lithium 123	1

Cat#	Description	Size	Batteries	Box
5-0237-91	Penlite Handle	143mm	2 "AA"	1
5-0237-92	Stubby Handle	110mm	2 "AA"	1
5-0237-93	Medium Handle	145mm	2 "C"	1
5-0237-96	Micro Handle	80mm	Lithium 123	1



Stainless Steel Howland Lock

- Natural lifting simplifies laryngeal exposure, making intubation possible even in the most difficult cases
- Built-in leverage decreases prying and reduces possibility of broken teeth
- Compatible with all hook-on conventional laryngoscope handles and blades
- Can often make an otherwise certain tracheotomy unnecessary

Cat#	Description	Box
5-0236-00	Conventional Howland Lock	1

Highest Quality. Unequaled Value. Outstanding Service.

(727) 548-4462 www.I-MA.com



Conventional Laryngoscopes Disposable LED Blades

SunMed's conventional stainless steel disposable blades answer the professional's request for a non-plastic disposable blade. These blades emit a clean, cool white LED light, and their low power use ensures a long battery life. The patented safety heel inhibits the blade tip from contacting handle. All sizes are compatible with Conventional ISO System handles. Using disposable blades helps reduce costs and concerns typically associated with reusable blades. These disposable blades are single use, sterile and individually packaged.

For detailed product information, visit www.Sun-Med.com or consult a SunMed product catalog.



Macintosh Standard / D™ LED

Miller Standard / D™ LED

Cat#	Size	Description	Box
5-5132-01	1	Infant	20
5-5132-02	2	Child	20
5-5132-03	3	Medium Adult	20
5-5132-04	4	Large Adult	20

Cat#	Size	Description	Box
5-5133-20	00	Premature	20
5-5133-00	0	Neonate	20
5-5133-01	1	Infant	20
5-5133-02	2	Child	20
5-5133-03	3	Medium Adult	20
5-5133-04	4	Large Adult	20

Highest Quality. Unequaled Value. Outstanding Service.

(727) 548-4462 www.I-MA.com