




Fiber Optic Laryngoscopes Greenline®/D™ Blades

SunMed's Fiber Optic Greenline® /D™ laryngoscope blades solve contamination problems and eliminate the cost and time spent cleaning blades and returning them for use. These stainless steel disposables also answer the professional's request for a non-plastic disposable blade. The polished acrylic stem produces exceptional illumination. All sizes are compatible with any reusable "Green System" Fiber Optic handle. These disposable blades are single use, sterile and individually packaged.

For detailed product information, visit www.Sun-Med.com or consult a SunMed product catalog.



Miller Greenline® / D™

Cat#	Size	Length	Description	Box
5-5333-20	00	67mm	Premature	20
5-5333-00	0	77mm	Neonate	20
5-5333-01	1	92mm	Infant	20
5-5333-02	2	100mm	Child	20
5-5333-03	3	130mm	Medium Adult	20
5-5333-04	4	155mm	Large Adult	20



Macintosh Greenline® / D™

Cat#	Size	Length	Description	Box
5-5332-01	1	92mm	Infant	20
5-5332-02	2	100mm	Child	20
5-5332-03	3	130mm	Medium Adult	20
5-5332-04	4	155mm	Large Adult	20



Robertshaw Greenline® / D™

Cat#	Size	Length	Description	Box
5-5337-01	0	94mm	Newborn	20
5-5337-02	0	112mm	Infant	20

Highest Quality. Unequaled Value. Outstanding Service.

(727) 548-4462 www.I-MA.com



Fiber Optic Laryngoscopes Handles

Choose from premium LED handles or standard xenon handles for use with SunMed's Fiber Optic laryngoscope blades. Both types are compatible with Green System blades.

For detailed product information, visit www.Sun-Med.com or consult a SunMed product catalog.



Penlite

Medium

Stubby



Large

Medium

Penlite

Stubby

Mini

Micro

Greenline SunBrite™ Fiber Optic LED Handle

- Complies with ISO 7376 Green Standard
- Eliminates any heat in handle
- Durable 75,000 hour LED lamp
- Cool, white illumination
- Batteries not included
- Low power use prolongs battery life

Standard Greenline® Handle

- Bright halogenated xenon lamp
- Economical
- Available in chrome plated brass and surgical stainless steel
- Compatible with all Green System Fiber Optic blades
- Batteries not included

Chrome Plated

Cat#	Description	Batteries	Box
5-0236-09	Medium Handle	2 "C"	1
5-0236-10	Penlite Handle	2 "AA"	1
5-0236-11	Stubby Handle	2 "AA"	1
5-0236-15	Large Handle	2 "D"	1

Surgical Stainless Steel

Cat#	Description	Batteries	Box
5-0236-39	Medium Handle	2 "C"	1
5-0236-30	Penlite Handle	2 "AA"	1
5-0236-31	Stubby Handle	2 "AA"	1
5-0236-35	Large Handle	2 "D"	1
5-0236-38	Mini Handle	Lithium "123"	1
5-0236-36	Micro Handle	Lithium "123"	1

Cat#	Description	Batteries	Box
5-0236-90	Penlite Handle	2 "AA"	1
5-0236-91	Medium Handle	2 "C"	1
5-0236-99	Stubby Handle	2 "AA"	1

Highest Quality. Unequaled Value. Outstanding Service.

(727) 548-4462 www.I-MA.com