

The Clear Choice for Disposable Infection Control



Cost-Effective
Lightweight
Comfortable

Timely Medical Innovations, LLC.

Since 1991 TMI has been providing disposable, personal protection products to help healthcare organizations comply with regulatory standards for infection control. Our focus on product innovation has yielded a line of low-cost, lightweight, comfortable eye and face protection products that have proven to be the preferred choice among healthcare professionals worldwide. As attention to personal protection and health safety grows, TMI will continue to be a market leader in low-cost, disposable personal protection products.

Resposables™

TMI's continuous focus on innovative product design has yielded many product patents over the years, the most recent of which is for our *Resposables™*. This line includes products for eye protection (*EyeShield™*) and face protection (*FaceShield™*), and has a unique, modular design consisting of reusable and disposable, optical grade lenses. *Resposables™* are lightweight, comfortable, latex-free and available in 10 designer colors, making them practical and stylish. The component frames and lenses are separately available to maximize cost-effective replenishment. And both products are CE certified for use in the European Union.



Resposables EyeShield™

- Extremely lightweight and comfortable (weighs only 3/10 of an ounce!)
- Patented design
- Latex-free
- "Wraparound" design provides both front and side eye protection
- Distortion-free, optical-grade lens
- Reusable frames available in 10 designer colors
- Fits comfortably with a surgical mask and over prescription eyewear
- UV-tinted and amber (blue-light absorbent) lenses also available
- Available in pre-assembled, ready-to-wear option
- Separately available replacement visors and lenses for maximized cost-effectiveness
- Meets the requirements for OSHA and EEC mandates



Resposables FaceShield™

- Extremely lightweight and comfortable (weighs less than 1.5 ounces!)
- Provides overhead, front and side eye protection
- 13" wide x 7" long distortion free, optical-grade shield with anti-static coating
- Latex-free
- Adjustable, and fog-resistant
- Reusable visors available in 10 designer colors
- Fits comfortably with a surgical mask and over prescription eyewear
- Includes a clear, snap-in nosepiece for support
- Amber (blue-light absorbent) shields also available
- Separately available replacement visors and lenses for maximized cost-effectiveness
- Meets the requirements for OSHA and EEC mandates



Resposables EyeShield Neck Straps™

Resposables EyeShield Neck Straps™ are the perfect accessory to drape the Eyeshield™ around one's neck to keep them in reach at all times .

- Reusable
- Available in 5 designer colors
- Latex-free



Disposable FaceShield™

The Disposable FaceShield™ is a one-piece disposable product for face protection.

- Standard Size (11.5" wide x 9" long)
- Full-Range Size (14" wide x 11.5" long) for added peripheral and extended face protection
- Distortion free, optical-grade lens with anti-static coating
- Lightweight, adjustable, and fog-resistant
- Latex-free
- Extra thick foam headband designed to comfortably conform to the forehead
- Fits comfortably with a surgical mask and over prescription eyewear



SideShield™

The reusable SideShield™ provides side-eye protection for wearers of prescription eyeglasses.

- Disposable clear lenses
- "Threads" easily onto existing eyewear
- Can be sterilized for re-use
- Conveniently packaged in reclosable bags



The Clip™

The Clip™ is a disposable nose clip used to seal airflow through the nose during pulmonary function testing and treatment.

- Comfortable non-slip, soft foam pads
- Snug fit
- Available in 10 designer colors