



Introducing the new KleanSpec® Disposable Vaginal Specula  
and Cordless Illumination System



An advancement in Women's Health  
technology that's nothing short of brilliant.

**WelchAllyn®**

Advancing Frontline Care™

# Now you can have better visualization of the exam area without the hassle of cords.

The new KleenSpec® Disposable Vaginal Speculum with the innovative Cordless Illumination System offers greater convenience for you, added comfort for your patients.


The Welch Allyn KleenSpec® Vaginal Specula and Illumination System has been totally redesigned to provide the key features you've come to expect—plus a number of important and exciting enhancements.

With so many new time-saving features and benefits, the new KleenSpec® Disposable Vaginal Specula and Cordless Illumination System are a must for every office and hospital.



The new 590 Series KleenSpec® Vaginal Speculum has all of the quality features you rely upon and more.

- > **Maneuverability**—new wider handle provides improved ergonomics, better balance, and it is much easier to manipulate during exams.
- > **Superior Illumination**—newly designed light pipe provides a homogeneous spot and eliminates glare back—all of the light is projected forward right where you need it.
- > **Comfortable**—smooth, molded acrylic will not pinch or bind tissue for greater patient comfort during the exam.
- > **Disposable**—Welch Allyn's single-use speculum eliminates the time and expense of reprocessing and reduces the chance of cross-contamination.
- > **Safe and Dependable**—stronger grip and the yoke is more flexible than low-cost substitutes.
- > **Compatibility**—the new 590 Series KleenSpec® Vaginal Speculum is compatible with the current illumination system (#78810) as well as the new illumination system (#79910).



#59000



#59001



#59004

590 Series KleenSpec® Vaginal Speculum  
Available in small (#59000), medium (#59001) and large (#59004) specula







The new Welch Allyn KleenSpec® Disposable Vaginal Specula 790 Series Illumination System will change the way you see your patients.

- > **Superior Visualization**—advanced LED light source provides better visualization of the exam area by supplying brighter, whiter light, right where it is needed.
- > **Cordless and Rechargeable**—no cord to get in the way or break during procedures, continuous on time of 80 minutes, recharge time of six hours, lithium ion battery technology. The new cordless illumination system easily fits right into the new KleenSpec® Disposable Vaginal Speculum.
- > **Automatic On/Off**—no more hassle with turning the unit on and off. When the new KleenSpec® Illumination System is inserted into the 590 Vaginal Speculum it is automatically illuminated, and once it is removed the light goes off.
- > **Versatility**—use your new Illumination System as an all-purpose light source as well as with the KleenSpec® Disposable Vaginal Speculum.
- > **Convenient**—no cords to clean, much less heat than traditional lamps.



**790 Series Illumination System**

Shown with #79920 KleenSpec® Disposable Vaginal Specula Cordless Illuminator with Charging Station and Auxiliary Power Cord

# The new KleenSpec® Disposable Vaginal Specula and Cordless Illumination System

## Ordering Information

### KleenSpec® Vaginal Specula—Premium 590 Series

- 59000 KleenSpec® 590 Series Premium Disposable Vaginal Specula, Small (24/box; 4 boxes/case)
- 59001 KleenSpec® 590 Series Premium Disposable Vaginal Specula, Medium (24/box; 4 boxes/case)
- 59004 KleenSpec® 590 Series Premium Disposable Vaginal Specula, Large (19/box; 4 boxes/case)

### KleenSpec® Box Series

- 59000-B KleenSpec® 590 Series Premium Disposable Vaginal Specula, Small (24/box)
- 59001-B KleenSpec® 590 Series Premium Disposable Vaginal Specula, Medium (24/box)
- 59004-B KleenSpec® 590 Series Premium Disposable Vaginal Specula, Large (19/box)

### KleenSpec® Vaginal Specula—Premium 590 Series w/Smoke Tube

- 59005 KleenSpec® 590 Series Premium Disposable Vaginal Specula with Smoke Tube, Small (12/box; 4 boxes/case)
- 59006 KleenSpec® 590 Series Premium Disposable Vaginal Specula with Smoke Tube, Medium (12/box; 4 boxes/case)

### KleenSpec® Disposable Vaginal Specula Illumination System

- 79910 Complete KleenSpec® Disposable Vaginal Specula Cordless Illuminator with Charging Station
- 79915 Same as above, for export purposes—specify country and voltage; conforms to IEC 601-1
- 79920 Complete KleenSpec® Disposable Vaginal Specula Cordless Illuminator with Charging Station and Auxiliary Power Cord
- 79925 Same as above, for export purposes—specify country and voltage; conforms to IEC 601-1
- 79900 KleenSpec® Disposable Vaginal Specula Cordless Illuminator
- 73910 Charging Station for Cordless Illuminator;  
For export purposes: 73912 (European), 73914 (UK), 73916 (Australia)
- 75910 Auxiliary Power Cord Accessory for Cordless Illuminator;  
For export purposes: 75912 (European), 75914 (UK), 75916 (Australia)

In addition, we offer a complete line of Vaginal Specula and Accessories to meet your needs:

### KleenSpec® Vaginal Specula—Premium 580 Series w/Sheath

- 58000S KleenSpec® 580 Series Premium Disposable Vaginal Specula with Sheath, Small (20/box; 4 boxes/case)
- 58001S KleenSpec® 580 Series Premium Disposable Vaginal Specula with Sheath, Medium (20/box; 4 boxes/case)
- 58004S KleenSpec® 580 Series Premium Disposable Vaginal Specula with Sheath, Large (15/box; 4 boxes/case)

### KleenSpec® Box Series

- 58000S-B KleenSpec® 580 Series Premium Disposable Vaginal Specula w/Sheath, Small (20/box)
- 58001S-B KleenSpec® 580 Series Premium Disposable Vaginal Specula w/Sheath, Medium (20/box)

### KleenSpec® Vaginal Specula—Economy Series

- 58600 KleenSpec® 586 Series Economy Disposable Vaginal Specula, Small (25/box; 4 boxes/case)
- 58601 KleenSpec® 586 Series Economy Disposable Vaginal Specula, Medium (25/box; 4 boxes/case)
- 78600 Expendable Illuminator (6/box; 6 boxes/case)—Can be sterilized 50 times before disposal

### KleenSpec® Disposable Vaginal Specula Illumination System

- 78810 KleenSpec® Illumination System—4.7 V Halogen Light Handle with on/off switch and transformer. 110-130 V, 60 Hz for North America. UL Listed, CSA Certified
- 78815 Same as above, for export purposes—Specify country and voltage; conforms to IEC 601-1
- 78800 Permanent 4.7 V Halogen Light Handle with On/Off Switch
- 73410 Transformer with 5' Cord—110-130 V, 50 Hz. UL Listed, CSA Certified
- 73315 Same as 73410, but for export purposes—Specify country and voltage; conforms to IEC 601-1
- 08800-U 4.6 V to 4.7 V Halogen Replacement Lamp for 78810 Illumination System
- 07800-U 6 V Halogen Replacement Lamp for 78010 Illumination System

### KleenSpec® Storage Systems

- 78820 KleenSpec® Illumination Storage System
- 58014 KleenSpec® Vaginal Speculum Dispenser

Patents Pending



Distributed By:

**Independent Medical Associates**

**IMA**  
Independent Medical Associates

11733 66th Street North  
Suite 113  
Largo, FL 33773

Phone: 888-548-4462 Fax: 888-548-1462  
E-Mail: Info@I-MA.com

**WelchAllyn®**

Advancing Frontline Care™