

Anesthesia Masks

Ultra Thin, Pliable Premium Contoured Masks

The
Clear
Choice



IMA
Independent Medical Associates

Anesthesia Masks

Ultra Thin, Pliable Premium Contoured Masks

All Westmed anesthesia masks feature clear, flexible faceplates and soft elastomeric cushions for maximum patient comfort.

- Eliminates the possibility of Corneal Abrasion due to contour of mask; Requiring less facial pressure.
- Unique Ergonomic Shape fits all facial configurations for improved sealing capabilities.
- Benefit to clinicians include; Requires less hand pressure, reduces strain on hands, and prevents hand discomfort.
- Soft, pliable air filled cushion with and without an inflation valve.
- Versatility of size — Includes six standard sizes to properly fit all your patients needs.
- The clear choice for patient and clinician comfort.



Part Number	Description	Case Quantity
-------------	-------------	---------------

700 Series - Ultra thin rotational molded cushion		
7-9700	Large adult w/hook ring and inflation valve	30
7-9701	Medium adult w/hook ring and inflation valve	30
7-9702	Small adult w/hook ring and inflation valve	20
7-9703	Pediatric w/hook ring and inflation valve	20
7-9704	Infant with inflation valve.	20
7-9705	Neonate with inflation valve	20

800 Series - Durable blow molded cushion		
7-9800	Large adult w/hook ring	30
7-9801	Medium adult w/hook ring.	30
7-9802	Small adult w/hook ring	20
7-9803	Pediatric w/hook ring	20
7-9804	Infant	20
7-9805	Neonate.	20

810 Series - Durable blow molded cushion with inflation valve		
7-9810	Large adult w/hook ring and inflation valve	30
7-9811	Medium adult w/hook ring and inflation valve. . . .	30
7-9812	Small adult w/hook ring and inflation valve	20
7-9813	Pediatric w/hook ring and inflation valve	20
7-9814	Infant with inflation valve.	20
7-9815	Neonate with inflation valve	20

900 Series - Durable blow molded BLUE cushion with extended back plate		
7-9900	Large adult w/hook ring and inflation valve	30
7-9901	Medium adult w/hook ring and inflation valve. . . .	30
7-9902	Small adult w/hook ring and inflation valve	20
7-9903	Pediatric w/hook ring and inflation valve	20
7-9904	Infant with inflation valve.	20
7-9905	Neonate with inflation valve	20

All masks are individually bagged.

Ask your local representative about other standard and custom configurations that are available to meet your anesthesia mask needs.



Independent Medical Associates
11733 66th St. N #113
Largo, FL 33771
Phone: 888.548.4462
Fax: 888.548.1462
www.I-MA.com