

Hyperinflation Sets

TOTAL CONTROL AT YOUR FINGERTIPS

Complete
Control



IMA
Independent Medical Associates

Hyperinflation Sets

TOTAL CONTROL AT YOUR FINGERTIPS

Westmed hyperinflation elbows and sets are precision engineered for versatility and functionality. They are the perfect choice for the NICU or post-operative transport.

- The high profile, knurled adjustment knob places full pressure control at the clinician's fingertips.
- A 180 degree swivel manometer port for superior mobility in confined spaces.
- ISO dimensions to ensure a secure, leak resistant attachment.
- Precision molded bleed orifice minimizes inadvertent PEEP or CPAP, while providing safe, effective pressure-limitation.



Hyperinflation Sets

Part Number	Description	Case Quantity
7-N7904-62	1/4 Latex Free bag, 84" O ₂ & pressure lines	12
7-N7905	1/2L Latex bag, 84" O ₂ & pressure lines	12
7-N7905-62	1/2 Latex Free bag, 84" O ₂ & pressure lines	12
7-N7901	1L Latex bag, 84" O ₂ & pressure lines	12
7-N7901-62	1 Latex Free bag, 84" O ₂ & pressure lines	12
7-N7902-62	2 Latex Free bag, 84" O ₂ & pressure lines	12
7-N7903-62	3 Latex Free bag, 84" O ₂ & pressure lines	12
7-N7910	Hyperinflation elbow, 84" O ₂ & pressure lines	50

Accessories

Part Number	Description	Case Quantity
0497	Manometer Cover	1
0499	Manometer	1

Ask your local representative about the no-charge Manometer program and other custom configurations that are available to meet your hyperinflation needs.

IMA
Independent Medical Associates

Independent Medical Associates
11733 66th St N #113
Largo, FL 33771
Phone: 888.548.4462
Fax: 888.548.1462
www.I-MA.com