

Oxygen Tubing  
Masks  
Water Trap



## Oxygen Tubing

- Kink resistant six channel inner lumendesign
- Available in two durometeroptions:  
75 shore (clear) and 85 shore (green)

75 durometertubingis soft, flexible tubingthat remains straightwhen uncoileddue to “no memory” design

85 durometertubinginhibitstwistingof the tubingwhich sometimesoccurs in longer lengthsdue to patientmovement



## Aerosol “Tee”

- Multi-use, universal respiratory connector(22 mmOD, 15 mmID)
- Loop on “T” facilitates hanging



## Super Spike Mask


It’s “cartoon” like appearance can give childrena calmingeffect when takingan aerosol respiratory treatment.

\*Sold with the VixOne™ and UniHeart™ nebulizer kits.



## Oxygen and Aerosol Masks

- Swivel snout accommodates multiple patient positions.
- Transparent mask allows clear view for patient assessment.
- Elastic strap with end caps are adjustableand easy to thread.



## In-Line Water Trap

- Helps eliminateunsafe and unnecessary moisture accumulation in the oxygen tubing
- Barbed connectors on end caps help prevent unwanted separation of tubing fromend caps



Product Number	Description	Case Quantity
<b>Oxygen Tubing</b>		
0007	7' Kink Resistant Tubing	50
0008	7' Kink Resistant Tubing w/threaded nut	50
0014	14' Kink Resistant Tubing	50
0022	21' Kink Resistant Tubing	25
0025	25' Kink Resistant Tubing	25
0027	25' Kink Resistant Tubing (Green)	25
0035	25' Kink Resistant Tubing (Purple)	25
0040	40' Kink Resistant Tubing	20
0050	50' Kink Resistant Tubing	20
0052	50' Kink Resistant Tubing (Green)	20
0055	50' Kink Resistant Tubing (Purple)	20
<b>Respiratory Accessories</b>		
0043	Oxygen Tubing Connectors (barbed)	50
0219	Aerosol Tee (Individually Bagged)	50
0221	In Line Water Trap for Oxygen Delivery	25
0407	Pediatric Tee	50
0409	Nipple and Nut (christmas tree)	50
0412	No-Smoking Sign	25
0480	Dry Humidifier, 6 psi	50
0497	Manometer Cover	1
0499	Manometer	1
0700	Air Cylinder Wrench	50
0704	Air Cylinder Wrench w/50 attachable chains	50
<b>Aerosol and Oxygen Masks</b>		
0250	Adult Aerosol Mask (Vented)	50
0251	Pediatric Aerosol Mask (Vented)	50
0360	Trach Mask, Adult	50
0375	Super Spike Pediatric Mask	50
0361	Trach Mask, Ped	50
0365	Adult Non-Rebreather Mask w/safety vent, 7' kink resistant tubing	50
0366	Pediatric Non-Rebreather Mask w/safety vent, 7' kink resistant tubing	50
0367	Adult Non-Rebreather Mask w/out safety vent, 7' kink resistant tubing	50
0370	Adult Medium Concentration Mask w/ 7' kink resistant tubing	50
0371	Pediatric Medium Concentration Mask w/ 7' kink resistant tubing	50
0760	Brass Regulator Washers	25