

New Product Update

Date: 06.24.09

To: Westmed Sales

From: Marketing

Super Spike Pediatric Aerosol Mask – Final Design

Our new Super Spike Mask has arrived! As part of the design change, we have placed the vent hole towards the bottom of the mask. This is to prevent medication from blowing directly into the patient's eyes. The freckles add character and give Super Spike more of a "puppy" look. All future shipments have the final design. As a reminder, I have listed the product codes Westmed offers with the Super Spike Mask.



Product Code	Description	Case Qty
0312	VixOne, 7' tubing, Super Spike Mask	50/cs
0313	VixOne, 7' tubing w/threaded nut, Super Spike Mask	50/cs
0375	Super Spike Mask	50/cs
100614A	Mini Heart Hi-Flo, 7' tubing, Super Spike Mask	25/cs
100628	Large Heart, 8' Corrugated Tubing, Super Spike Mask	18/cs
100942	UniHeart, 7' tubing, Super Spike Mask	25/cs

IMA
Independent Medical Associates