


Grip-Lok® H is available in two sizes to meet most pediatric and neonatal applications.

Shown Actual Size



2.5" (6.4 cm)



3" (7.6 cm)


Product Number: **3100ESH**
Size: Extra Small
Tubing Range: 1.0-3.0mm

Product Number: **3200SH**
Size: Small
Tubing Range: 1.5-4.5mm



New Grip-Lok H release liners make it easy to apply.

Grip-Lok® NGH is ideal for feeding tube securement applications.



Grip-Lok NGH is easy to apply even with gloved hands.



1.625" (4.1 cm)

Product Number: **2100NGH**
Size: Small
(can be cut in half prior to removal of liners to reduce size)
Tubing Range: 1.0-2.0mm



**Safe
Simple
Secure
Gentle**



New!
Grip-Lok NGH
for Foley Tube
Catheters



IMA

Independent Medical Associates

11733 66th St N #113; Largo, FL 33773
P: 888-548-4462; F: 888-548-1462
info@I-MA.com; www.I-MA.com

GRIP-LOK® H

*Hydrocolloid Tubing
Securement Device Designed
Specially for Pediatric
& Neonatal Patients*